

Feeder Seat System Sizing Chart				
Dimension	Small Feeder Seat	Medium Feeder Seat	Large Feeder Seat	XLarge Feeder Seat
Overall Height	21" (53 cm)	27" (68.5 cm)	34" (86 cm)	41" (104 cm)
Outside Width	10 ¹ / ₂ " (26.5 cm)	14" (35.5 cm)	17" (43 cm)	21" (53 cm)
Inside Depth	8" (20 cm)	9 ¹ / ₂ " (24 cm)	13 ¹ / ₄ " (33 cm)	16 ¹ / ₂ " (42 cm)
Inner Width	7 ³ / ₄ " (18.7 cm)	9" (22.8 cm)	12 ¹ / ₂ " (31.7 cm)	15" (38 cm)
Inner Height	18" (46 cm)	23" (58.5 cm)	29" (73 cm)	35" (89 cm)
Child's Height	up to 36" (91 cm)	up to 48" (122 cm)	up to 60" (150 cm)	up to 72" (183 cm)

The Feeder Seat System provides secure positioning and comfort during meals and activity.