



DropSafe® Safety Pen Needles

Intelligent design with viewing window for easy priming, while improving patient comfort.1



You can maintain a high standard of safety, efficacy, and comfort for clinicians and patients with DropSafe®.

- Improve patient comfort with the first 31G safety pen needle available in North America, DropSafe 31G.*
- In a simulated clinical study evaluating the safety and performance of DropSafe in 540 end-to-end patient usage trials, no failures were observed and no needle stick injuries occurred.²
- Prevent accidental needle sticks with DropSafe's patient end needle shield. It is recommended to cover back-end with outer cap & dispose. Put the outer cap on the needle back-end.3
- See the needle tip clearly during priming with DropSafe's needle viewing window.¹







Needle Viewing Window

Clear, unobstructed visibility to needle tip during priming allows caregiver to see the droplet of primed drug.¹



Thoughtful Safety Design

A needle shield prevents accidental needle sticks at the patient end. It is recommended to cover back-end with outer cap & dispose. Put the outer cap on the needle back-end.³



Unique In-house Lubrication

Droplicon™ needle lubrication designed for smooth injection.



Needle Lock-out Confirmation

Needle is automatically locked after use to prevent re-use. Red stripes appear to confirm needle has been used and is ready to be disposed of safely.²



* As of December 31, 2018. **1**. HTL-STREFA S.A.

Data on file. Qualitative study on safety pen needles – Focused group interviews with nurses in Spain, PMR.

2018. 2. HTL-STREFA S.A. Data on file. Final report

Simulated clinical use testing on the DropSafe Safety
 Pen Needle. Version 1.0. 2017.
 3. HTL-STREFA S.A.

Data on file. Instructions For Use - DropSafe Safety

Ordering Information:

| MFR# | NDC# | Pr | Product Size | | | Color | Quantities | |
|------|--------------|----|----------------------------|------------|--------|--------|-----------------------|--|
| 8162 | 8489-8162-10 | 31 | G (0.2 | 5mm) × 6mn | n (¼") | Purple | 100 pieces per box | |
| 8163 | 8489-8163-10 | 31 | 31G (0.25mm) × 8mm (5/16") | | | Blue | 2,000 pieces per case | |

