

# ETC Line, Five Drawer, Treatment Cart #ETC-5



### Over 20 Available Colors:

Click or copy & paste the link below to view all current available color options:

https://www.harloff.com/finish-options/

650 Ford Street Colorado Springs, CO 80915 Toll:(800)433-4064 \* Ph: (719)637-0300 Fax:(719)597-8273 \* www.harloff.com

# **Specifications**

#### **Construction:**

• Sturdy, uni-body steel construction

#### Features:

- Low height make these carts ideal for bedside or under a counter
- · Key lock with two keys
- 3" stem casters, two with brake
- · Roller drawer slides
- · Removable plastic top
- Latex free
- Twelve year warranty
- · Hard-baked powder coat finish

## **Drawer Configuration:**

- (3) 2.75" (7.0 cm)
- (2) 5.75" (14.6 cm)

#### **Assembled Cart Dimensions**

- 28.63" H x 23.38" W x 15.13" D (72.73 cm H x 59.4 W x 38.4 D)
- Cart weight: 75 lbs (34.1 kg)
- Approx. Shipping Weight: 81 pounds (36.7 kg)