

ENCORE® Latex Textured

SURGICAL GLOVES

Latex, powder-free surgical gloves ensuring secure grip for confident instrument handling

- Enhanced grip: The textured finish of these sterile latex gloves delivers exceptional grip of medical instruments under dry or damp conditions, making them ideal outer gloves when double gloving
- Heightened tactile sensitivity: The unique latex formulation and processes enable the thickness of these sterile gloves to be optimized for the wearer to experience the ideal balance of tactile sensitivity and hand protection
- Effortless donning and doffing: ENCORE[®] Latex Textured surgical gloves are designed with a proprietary internal coating to facilitate easy dry and damp donning when first worn or during intraoperative changes
- Chemo-tested: Tested to ASTM D6978, these textured latex gloves are safe to use with a wide range of chemotherapy drugs*

CAUTION: Safe use of these gloves by latex-sensitized individuals has not been established. This product contains natural rubber latex, which may cause allergic reactions.

*Contact Ansell Customer Service for specific chemotherapy drug permeation times and recommendations

Recommended For

- Enhanced grip in wet environments
- Extra comfort during long, sustained procedures
- Damp donning
- Use as an outer glove when double gloving

FEATURED TECHNOLOGIES











ENCORE® The Confidence to Perform

Dedicated to delivering optimized solutions engineered to provide a high degree of confidence to surgeons, O.R. nurses and surgical technicians seeking protection with a high level of dexterity, comfort and grip so they can perform at their best.



ENCORE® Latex Textured

SURGICAL GLOVES

Material	Natural Rubber Latex							
Powder Content	Powder-Free							
Color	Natural							
Thickness	Standard							
Cuff Style	Beaded with SUREFIT™ Technology							
External Glove Surface	Textured Fingers and Palm							
Internal Glove Surface	Polymer Coated with DERMASHIELD™ Technology							
Allergy Prevention	None							
Grip Level	High							
Double Gloving Recommendation	Outer glove							
Tested For Use With Chemotherapy Drugs	Yes, in accordance with ASTM D6978 (Not listed in the US FDA 510k)							
Freedom from Holes (Inspection level I)	0.65 AQL							
Sterile	Yes							
Product Standards	ASTM D3577, ASTM D6978, ASTM D7160, ASTM F739, EN 16523-1, EN 374-2 and -4, EN 455 1-4 EN ISO 21420:2020, EN ISO 374-1:2016, EN ISO 374-5:2016, ISO 10282							
Quality / Environmental Standards	EN 556, ISO 11137-Part 1, ISO 13485, ISO 14001, ISO 9001							
Packaging	50 pairs per box; 4 boxes per carton/case; 200 pairs per carton/case							
Shelf Life	3 years							
Storage Instructions	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.							
Vulcanization Chemical Accelerators	Zinc Diethy	rldithiocarb	amate (ZDEC	C), Zinc Dibu	tyldithiocark	oamate (ZDE	BC)	
ORDERING INFORMATIO	N							
Size	5.5	6	6.5	7	7.5	8	8.5	9
Product Code	5785000	5785001	5785002	5785003	5785004	5785005	5785006	578500

Contact your Ansell representative for ordering or more information.

North America US Tel: 800 952 9916 CA Tel: 844 494 7854 Email: insidesalesus@ansell.com Sandel US Tel: 866 764 3327 Sandel email: sandel-cs@ansell.com Central & South America Tel: +52(442) 296 20 50 Email: cslac@ansell.com Brazil Tel: +55-11-3356-3100 Email: luvas.medicas@ansell.com Europe, Middle East & Africa Tel: +32 (0) 2 528 74 00 Email: info@ansell.eu

Asia Pacific Tel: +603 8310 6688 Email: apac.medical@ansell.com India
Tel: +91 22 38124500
Email: info.india@ansell.com
Australia & New Zealand
Tel: +61 3 9270 7270
Email: protection@ap.ansell.com
China
Tel: +86(21) 3827 5005

Email: infochina@ansell.com