


iontophoresis systems and pain relieving patches

IontoPatch®

IontoPatch™



IontoPatch SP™



	IontoPatch® 80	IontoPatch® SP	IontoPatch® STAT	IontoPatch Extra Strength™
active area	15.5 cm ²	8.0 cm ²	12.0 cm ²	12.0 cm ²
dosage	80 mA-min	40 mA-min	80 mA-min	120 mA-min
fill volume	1.2 - 1.6 cc	1.2 - 1.6 cc	1.2 - 1.4 cc	1.2 - 1.6 cc
wear time	14 hours	14 hours	4 hours	8 hours

- self-contained wireless technology eliminates need for external batteries or wires
- allows patient to return to daily activities while receiving time-released iontophoresis
- made of flexible, breathable material
- technician applies medication directly to the electrode
- shuts off automatically when prescribed dose has been administered
- 6 treatment kits per box (patch, saline vial, alcohol prep)

	description	qty	
13-5220	IontoPatch® 80	16 ea	60.00
13-5221	IontoPatch® SP	16 ea	60.00
13-5223	IontoPatch® Stat	16 ea	75.00
13-5224	IontoPatch Extra Strength™	16 ea	85.00

Iomed Companion 80™

- reserve battery capacity: compensates for patients with higher skin resistance
- hypoallergenic adhesive: adheres well, even in the shower
- high-tech design conforms well to treatment areas and fits comfortably beneath clothing
- 6 treatment kits per box (patch, saline vial, alcohol prep)



Companion 80™	
active area	6.45 cm ²
dosage	80 mA-min
fill volume	1.1 cc
wear time	24 hours





	description	qty	
50-0004-4	Iomed Companion 80™	6 ea	52.50

treatment modalities

LidoFlex™ pain relieving patches

- pain relief patch that delivers clinical strength lidocaine with the highest percentage available without a prescription (4%)
- topical anesthetic/analgesic patch targets specific pain with deliberate, localized relief
- unique form-fitting shapes can be cut to size, are easily applied/removed and are waterproof
- non-habit forming
- available in 5-patch or 20-patch packs



pack with 5 patches each	case quantity 30 packs with 5 patches each	MAP <i>(minimum advertised price applies)</i>	pack with 20 patches each	case quantity 15 packs with 20 patches each
13-5410-5 32.50	13-5410-5-30 975.00	knee	13-5410-20 117.50	13-5410-20-15 1,762.50
13-5411-5 32.50	13-5411-5-30 975.00	elbow	13-5411-20 117.50	13-5411-20-15 1,762.50
13-5412-5 32.50	13-5412-5-30 975.00	back	13-5412-20 117.50	13-5412-20-15 1,762.50
13-5413-5 32.50	13-5413-5-30 975.00	shoulder	13-5413-20 117.50	13-5413-20-15 1,762.50
13-5414-5 32.50	13-5414-5-30 975.00	heel	13-5414-20 117.50	13-5414-20-15 1,762.50
13-5415-5 32.50	13-5415-5-30 975.00	flex strip, 2/ea	13-5415-20 117.50	13-5415-20-15 1,762.50