

- the Baseline[®] hydraulic hand dynamometer gives accurate grip strength readings without the subject being able to “feel” the handle move
- readings are repeatable
- maximum reading remains until reset
- readings in lb and kg
- CE certified
- manufactured and serviced in the USA since 1978

- recalibration available once/year based on warranty
- five position handle and body contours assure results consistent with published Baseline[®] and Jamar[®] studies
- internationally accepted design assures reliability, accuracy and repeatability
- lightweight units include protective carrying case



200 lb capacity

- inexpensive and reliable
- 200 lb (90 kg) capacity
- 2½" diameter dial
- attractive red finish



LITE[®]
1 year warranty

12-0241 200 lb LITE[®] 225.00

- our most popular style
- 200 lb (90 kg) capacity
- 2½" diameter dial
- blue hammertone finish



Standard
2 year warranty

12-0240 200 lb Standard 325.00

- maximum warranty
- 200 lb (90 kg) capacity
- 2½" diameter dial
- black finish



HD[®]
5 year warranty

12-0221 200 lb HD[®] 425.00

extended range 300 lb capacity

- extra large 3½" diameter dial
- 200 lb (90 kg) capacity
- blue hammertone finish



Hi-Res[®]
2 year warranty

12-0243 200 lb HiRes[®] 350.00

- extended range 300 lb (136 kg) capacity
- extra large 3½" diameter dial
- blue hammertone finish



Hi-Res[®]
ER[™]
2 year warranty

12-0246 300 lb ER[™] HiRes[®] 375.00



- extended range 300 lb (136 kg) capacity
- electronic zero calibration system
- max clear button
- toggle maximum and current reading
- toggle lb / kg
- 2½" diameter LCD gauge **digital**
- 2 AAA batteries included **2 year warranty**
- blue hammertone finish

12-0247 300 lb digital 725.00