

Baseline[®] spirometer



spirometer with plastic mouth pieces

- compact, lightweight spirometer for screening large populations
- subject blows into disposable mouthpiece; stays at max until reset
- indicator registers forced vital lung capacity up to 7 liters



additional disposable mouthpieces


spirometer

12-1710	spirometer w/ 50 plastic mouthpieces	150.00
---------	--------------------------------------	--------

spirometer accessories

12-1711	100 plastic mouthpieces	25.00
12-1714	250 disposable mouthpieces	25.00
12-1712	500 disposable mouthpieces	50.00
12-1713	1000 disposable mouthpieces	75.00

Baseline[®] body fat scale



- tempered glass scale measures weight, body fat %, water %, muscle %, bone density
- features an LCD screen, auto on, 5 activity levels and memory for 12 users
- 330 lb/150 kg capacity
- 4% - 45% body fat; 37.8% - 66% body water
- 0.2 lb/0.1% BF/water resolution
- 12-person memory
- 12½"L x 12¼"W x 2"H
- 3V battery included

12-1190	body fat scale	45.00
---------	----------------	-------

Omron[®] body monitor scales

HBF-510W body composition

- measures body fat %, BMI, weight, skeletal muscle and visceral fat
- 4-person memory
- 4 AA batteries included



12-1193	Omron [®] HBF-510W	80.00
---------	-----------------------------	-------

HBF-400 fat loss monitor

- displays weight, BMI and body fat %
- 4-person memory
- 4 AA batteries included




12-1192	Omron [®] HBF-400	60.00
---------	----------------------------	-------

Baseline[®] Sit 'n Reach[®] trunk flexibility test

- perform sit 'n reach flexibility tests
- measurements include: fitness, program, Cooper Institute/YMCA, SHAPE America (formally known as AAHPERD)
- scales printed in cm and inches
- sturdy metal construction with a built-in foot plate and a maximum stretch indicator that remains until reset
- deluxe model offers adjustable start point
- deluxe model can be used to normalize scores to account for variation in the arm and leg length of the person being tested



standard model



deluxe model

with adjustable start point measurement mechanism

12-1085	standard	125.00
12-1086	deluxe	200.00

Fabrication Enterprises, PO Box 1500, White Plains, New York 10602 USA

tel: 800-431-2830 • 914-345-9300 Orders@FabEnt.com fax: 800-634-5370 • 914-345-9800