

ivNow®



Customize to fit your warming needs.

The ivNow can be ordered in a 1, 2 or 3 pod configuration.

HCAI - OPM-0060-13



- Know the warming time and temperature of every IV bag.
- Warms 1 liter iv bags in 30 minutes or less.
- Point-of-care use saves space, time & footsteps.
- Requires no disposables.



ivNow-1

- Capacity: 1 pod (0.5 to 3 liter bags)
- Exterior Size (HxWxD): 7.1" (180mm) x 13.2" (335mm) x 8" (203mm)



ivNow-2

- Capacity: 2 pods (1 to 6 liter bags)
- Exterior Size (HxWxD): 14.2" (360mm) x 13.2" (335mm) x 8" (203mm)



ivNow-3

- Capacity: 3 pods (1.5 to 9 liter bags)
- Exterior Size (HxWxD): 21.2" (540mm) x 13.2" (335mm) x 8" (203mm)



ivNext®

Inventory Management System



Intelligent & Innovative Fluid Warming

- Know when an IV bag has reached the recommended set point temperature (adjustable injection mode)
- Color-coded lighting indicates bag readiness or expiration.
- Organized inventory allows for faster retrieval of fluids.

*Check with your local distributor for availability if purchasing outside of the US.



- Total capacity: 10 1-liter bags
- Interior size: 2.5 ft³
- Exterior Size (HxWxD): 22" (559mm) x 18.5" (470mm) x 23.5" (597mm)
- Stackable

STORAGE CABINET



1560-GS

- Interior size: 15.6 ft³
- Exterior Size (HxWxD): 75" (1905mm) x 30" (762mm) x 21.25" (540mm)

Pair with a Combination Warmer
Shown with EC1350BL



MAKE IT MOBILE

MS-400 Cart

- Mobile Cart for small blanket and fluid warmers
- 36.5" (927mm) H; 26.0" (660mm) W; 26.0" (660mm) D
- Fits the following models: DC250/L, DC350, DC400/L, EC250/L, EC350/L, EC400/L, ivNext®



Shown with EC250L Fluid Warmer

MS-550 Cart

- Mobile Cart for multipurpose warmers
- 34" (866mm) H; 32.5" (825mm) W; 29.5" (749mm) D
- Fits the EC550BL-18" EC550-24" depth



Shown with EC550BL-18" Depth Combination Warmer



enthermics®
medical systems



INNOVATIVE & AFFORDABLE PATIENT WARMING SOLUTIONS



5-Year Parts, 1-Year Labor Warranty on cabinets
*U.S./Canada only

Entermics is committed to providing the best value and return on investment for our customers. Our never ending drive for innovation and quality has made Entermics a leader in patient warming. As an ISO 13485:2016 certified company, we provide sales and technical support to our customers worldwide. Each warming unit is proudly manufactured and inspected in the U.S.A.



ENTHERMICS MEDICAL SYSTEMS

An ISO 13485:2016 Certified Company

PO Box 443, Menomonee Falls WI 53052-0443

W164 N9221 Water St, Menomonee Falls, WI 53051

Phone 262-251-8356 | 800-862-9276

generalinfo@enthermics.com | www.enthermics.com

Made in the USA

TITAN SERIES

Meet all your blanket & fluid warming needs with the Entermics Titan and Comfort Series Warmers.

WarmSafe™ zone heat – warms safely, quietly & evenly by adjusting temperature throughout the cabinet

Trust the durable and affordable stainless steel, Titan Series warmers for efficient blanket & fluid warming.

Enjoy the convenient, patient-friendly Comfort Series warmers for point-of-care placement.

*Patent No: US 8,217,316; US 8,581,152

EC250

- Capacity: 4-6 blankets
- Interior size: 2.5 ft³
- Exterior size (HxWxD): 22" (559mm) x 18.5" (470mm) x 23.5" (597mm)
- Countertop model
- Stackable
- No fans or moving parts



EC350

- Capacity: 10-12 blankets
- Interior size: 3.5 ft³
- Exterior size (HxWxD): 28" (711mm) x 18.5" (470mm) x 23.5" (597mm)
- Floor or countertop model
- Stackable
- No fans or moving parts



EC400

- Capacity: 12-16 blankets
- Interior size: 4.0 ft³
- Exterior size (HxWxD): 22" (559mm) x 24.0" (610mm) x 25.75" (654mm)
- Stackable
- No fans or moving parts



DC250

- Capacity: 4-6 blankets
- Interior size: 2.5 ft³
- Exterior size (HxWxD): 22" (559mm) x 18.5" (470mm) x 22.0" (559mm)
- Countertop model
- Stackable
- No fans or moving parts



DC350

- Capacity: 10-12 blankets
- Interior size: 3.5 ft³
- Exterior size (HxWxD): 28" (711mm) x 18.5" (470mm) x 22.0" (559mm)
- Stackable
- No fans or moving parts



EC750

- Capacity: 20-25 blankets
- Interior size: 7.5 ft³
- Exterior size (HxWxD): 35" (889mm) x 24.0" (610mm) x 25.75" (654mm)
- Stackable
- No fans or moving parts



EC1540

- Capacity: 40-50 blankets or 60 1-liter irrigation bottles
- Interior size: 15.4 ft³, 7.7 ft³ per chamber
- Exterior size (HxWxD): 72.75" (1848mm) x 24" (610mm) x 32.5" (870mm)
- No fans or moving parts
- Shown with optional bumper



*Check with your local distributor for availability if purchasing outside of the US.

EC2060


- Capacity: 60-70 blankets
- Interior size: 20.6 ft³
- Exterior size (HxWxD): 71" (1803mm) x 32" (813mm) x 29.25" (743mm)
- HCAI - OPM-0067-13



*Check with your local distributor for availability if purchasing outside of the US.

DC400

- Capacity: 12-16 blankets
- Interior size: 4.0 ft³
- Exterior size (HxWxD): 22" (559mm) x 24" (610mm) x 24.3" (617mm)
- Stackable
- No fans or moving parts



DC750

- Capacity: 20-25 blankets
- Interior size: 7.5 ft³
- Exterior size (HxWxD): 35" (889mm) x 24.0" (610mm) x 24.3" (617mm)
- Stackable
- No fans or moving parts



BLANKET WARMERS

Provide patients with warm blankets and cabinets that are engineered for years of reliable operation.

EC250L

- Total capacity: 12 1-liter bottles or 16 1-liter bags
- Interior size: 2.5 ft³
- Exterior size (HxWxD): 22" (559mm) x 18.5" (470mm) x 23.5" (597mm)
- Counter top model
- Stackable



EC350L

- Total capacity: 24 1-liter bottles or bags
- Interior size: 3.5 ft³
- Exterior size (HxWxD): 28" (711mm) x 18.5" (470mm) x 23.5" (597mm)
- Floor or counter top model
- Stackable



EC400L

- Total capacity: 20 1-liter bottles or 28 1-liter bags
- Interior size: 4.0 ft³
- Exterior size (HxWxD): 22" (559mm) x 24.0" (610mm) x 25.75" (654mm)
- Floor or counter top model
- Stackable



DC250L

- Total capacity: 12 1-liter bottles or 16 1-liter bags
- Interior size: 2.5 ft³
- Exterior size (HxWxD): 22" (559mm) x 18.5" (470mm) x 22.0" (559mm)
- Counter top model
- Stackable



DC400L

- Total capacity: 20 1-liter bottles or 28 1-liter bags
- Interior size: 4.0 ft³
- Exterior size (HxWxD): 22" (559mm) x 24.0" (610mm) x 24.3" (617mm)
- Floor or counter top model
- Stackable



FLUID WARMERS

Meet healthcare guidelines and warm IV and irrigation fluids to recommended temperatures with WarmSafe™ heating technology.

MULTI-PURPOSE WARMERS

EC550BL-18

- Total capacity:
 - TriMode chamber: 24 1-liter bottles/ 20 1-liter bags or 14-18 blankets
- Exterior size (HxWxD): 25.12" (638mm) x 30" (762mm) x 21.12" (536mm)
- 18" chamber depth
- Also available as a built-in cabinet

Unit shown with fluid bottles, bags and blankets for illustration purposes only



EC550BL-24

- Total capacity:
 - TriMode chamber: 30 1-liter bottles + 36 0.5-liter bottles/ 48 1-liter bags or 16-20 blankets
- Exterior size (HxWxD): 25.12" (638mm) x 30" (762mm) x 28.5" (724mm)
- 24" chamber depth
- Also available as a built-in cabinet

Unit shown with optional basket - set up for fluids



COMBINATION WARMERS

EC1350BL

- Total capacity:
 - Blanket chamber: 24-30 blankets; 9.39 ft³ interior size
 - TriMode chamber: 24 1-liter bottles/ 20 1-liter bags or 9-12 blankets
- Interior size: 13.3 ft³
- Exterior size (HxWxD): 75" (1905mm) x 30" (762mm) x 21.25" (540mm)
- Available as a built-in or pass-through cabinet
- HCAI - OPM-0562-19



EC1850BL

- Total capacity:
 - Blanket chamber: 30-40 blankets; 13 ft³ interior size
 - TriMode chamber: 30 1-liter bottles/ 24 1-liter bags or 10-14 blankets
- Interior size: 18.5 ft³
- Exterior size (HxWxD): 75" (1905mm) x 30" (762mm) x 28.5" (724mm)
- HCAI - OPM-0067-13
- *Height with casters

