



MYCO

MEDICAL

TheSharpChoice®

NMSDC Certified Diversity Supplier

DATASHEET




Pictured: Glassvan® Stainless Steel Surgical Blade, Size #15

GLASSVAN® Surgical Blades

Glassvan® Surgical Blades have been trusted by users around the world for more than 25 years.

- Automated grinding and proprietary polishing processes ensure consistent, clean, and sharp cutting edges
- Produced from high-quality Swiss carbon and stainless steel

- Two VCI strips ensure product and packaging integrity and reduce risk of blade piercing the foil
- Manufactured in FDA-inspected facilities and conform to ISO-7740 and ISO-13485
- Carbon steel and stainless steel blades available sterile in sizes 10-25, with select sizes also available non-sterile
- Carbon steel miniature blades available in sizes 61, 62, 64, 67, 67MIS
- Podiatry blades available in stainless steel sizes 9, 17, 81, 84, 88, and in carbon steel sizes 312, 313, 314, 316




Sizes 10-25 available in Swiss carbon and stainless steels

BRAND (Company)	General Surgery Blades in Carbon and Stainless Steels	VCI Strips	Miniature & Podiatry Blades
Glassvan® Surgical Blades (MYCO Medical)	✓	✓	✓
Bard-Parker™ Surgical Blades (Aspen)	✓	✓	✓
Surgical Blades (Miltex)	✓		
Surgical Blades (Southmedic)	Stainless Steel Only		

Glassvan® Stainless Steel Surgical Blades

BLADE SIZE	STERILE PACKAGING		NON-STERILE PACKAGING	
	MYCO Code	Unit of Measure	MYCO Code	Unit of Measure
#10	3001T-10	100/bx, 10 bx/cs	3000-10*	500/bx, 20 bx/cs
#11	3001T-11	100/bx, 10 bx/cs	3000-11*	500/bx, 20 bx/cs
#12	3001T-12	100/bx, 10 bx/cs	3000-12*	500/bx, 20 bx/cs
#12B	3001T-12B	100/bx, 10 bx/cs	3000-12B*	500/bx, 20 bx/cs
#15	3001T-15	100/bx, 10 bx/cs	3000-15*	500/bx, 20 bx/cs
#15C	3001T-15C	100/bx, 10 bx/cs	3000-15C*	500/bx, 20 bx/cs
#20	3001T-20	100/bx, 10 bx/cs	3000-20*	500/bx, 20 bx/cs
#21	3001T-21	100/bx, 10 bx/cs	3000-21*	500/bx, 20 bx/cs
#22	3001T-22	100/bx, 10 bx/cs	3000-22*	500/bx, 20 bx/cs
#23	3001T-23	100/bx, 10 bx/cs	3000-23*	500/bx, 20 bx/cs
#24	3001T-24*	100/bx, 10 bx/cs	3000-24*	500/bx, 20 bx/cs
#25	3001T-25	100/bx, 10 bx/cs	3000-25*	500/bx, 20 bx/cs
#60	3001T-60	100/bx, 10 bx/cs	3000-60*	500/bx, 20 bx/cs

*Special Order Item — Minimum orders and lead times may apply



Sterile Glassvan® Carbon Steel Miniature Blades

BLADE SIZE	MYCO CODE	UNIT OF MEASURE
#6100	2002-61	12/bx, 50 bx/cs
#6200	2002-62	12/bx, 50 bx/cs
#6400	2002-64	12/bx, 50 bx/cs
#6700	2002-67	12/bx, 50 bx/cs
#67MIS	2002-67MIS	12/bx, 50 bx/cs



Consistent, clean edges

Sterile Glassvan® Carbon Steel Surgical Blades

BLADE SIZE	MYCO CODE	UNIT OF MEASURE
#10	2001T-10	100/bx, 10 bx/cs
#10A	2001T-10A*	100/bx, 10 bx/cs
#11	2001T-11	100/bx, 10 bx/cs
#12	2001T-12	100/bx, 10 bx/cs
#12B	2001T-12B	100/bx, 10 bx/cs
#15	2001T-15	100/bx, 10 bx/cs
#15C	2001T-15C	100/bx, 10 bx/cs
#20	2001T-20	100/bx, 10 bx/cs
#21	2001T-21	100/bx, 10 bx/cs
#22	2001T-22	100/bx, 10 bx/cs
#23	2001T-23	100/bx, 10 bx/cs
#24	2001T-24*	100/bx, 10 bx/cs
#25	2001T-25	100/bx, 10 bx/cs
#70	2001T-70GV	100/bx, 10 bx/cs
Stitch Cutter	2002-SC	100/bx

*Special Order Item — Minimum orders and lead times may apply

Sterile Glassvan® Podiatry Blades

BLADE SIZE	MATERIAL	MYCO CODE	UNIT OF MEASURE
#9	Stainless Steel	2003-09	12/bx, 50 bx/cs
#17	Stainless Steel	2003-17	12/bx, 50 bx/cs
#81	Stainless Steel	2003-81	12/bx, 50 bx/cs
#84	Stainless Steel	2003-84	12/bx, 50 bx/cs
#88	Stainless Steel	2003-88	12/bx, 50 bx/cs
#312	Carbon Steel	2003-312	12/bx, 50 bx/cs
#313	Carbon Steel	2003-313	12/bx, 50 bx/cs
#314	Carbon Steel	2003-314	12/bx, 50 bx/cs
#316	Carbon Steel	2003-316	12/bx, 50 bx/cs

As a leading NMSDC Diversity Certified Supplier, MYCO Medical is committed to “Bringing Reliability, Innovation and Value” to the Healthcare industry with an emphasis on “physician preference” surgical devices and disposables.

To order, learn more or receive complimentary evaluation samples:

Please call MYCO Medical at (800) 454-6926 or (919) 460-2535 or email us at info@mycomedical.com

See the complete line of MYCO innovative disposable and reusable medical devices at www.mycomedical.com



2015 Production Drive
Apex, NC 27539 USA

info@mycomedical.com
www.mycomedical.com

Some medical devices are not licensed for sale in Canada.

© 2017 MYCO Medical, Inc. All rights reserved.

LIT-00226 Rev 1 Eff 09.18.2017