## **INVACARE®**

# 9805 AND 9805P HYDRAULIC LIFTS





Too often, manual patient transfers lead to caregiver back injuries. The **Invacare® 9805 and 9805P Hydraulic Lifts** were created to help make handling situations safer and affordable for everyone involved. With their slim design, these hydraulic lifts are versatile enough to use for most any and all patient transfers. The features of the 9805 and 9805P include a padded swivel bar and push handle, a 450 lb. weight capacity, durable chrome-plated construction or painted finish and 360° patient rotation without side-to-side sway. Lightweight construction and easy disassembly allow for quick transport and setup. Indispensable to patients and caregivers alike, the 9805 and 9805P lifts offer safety, comfort and stability.

#### **Features**

Six-point swivel bar attachment easily adapts to all sling styles and positions

Low friction casters significantly improve rollability and contribute to caregiver safety

20"-64" range allows for patient to be picked up from lying position on floor

Pump handle can rotate from side to side for convenience of attendant

Offset mast and boom style provide better lift path, maximizing full travel range

Swivel bar and pump handle are padded for comfort Elongated handle design serves a broad range of operator heights and offers improved handling ability Completely enclosed cam moves legs of adjustable base Base offers obstruction-free operation; ideal for cleaning Screw-in pin at mast base provides maximum stability and

prevents removal of mast


## **HCPCS Code: E0630**











- 2. Hydraulic pump
- 3. Swivel bar
- 4. Standard sling

	-
Model #	Chain or Strap Assembly Kit
(required	for Invacare 9000 Series Wheelchairs)
9070	sling strap assembly kit
9071	sling chain assembly kit

Model #	Fabric/Style
9042	solid polyester without commode
9043	solid polyester/with commode
9046	mesh polyester/without commode
9047	mesh polyester/with commode







Base



Casters



Weight





Warranty

Model # 9805, 9805P

At Sling Hook Up Dimensions 20"-64" H

6.5"H Clearance 42.25"W Open 22"W Closed 46.75"L

5"

68 lb.

Weight Capacity 450 lb.

3 Years Lift 1 Year Hydraulic Pump

# **Invacare Corporation**

www.invacare.com

### **USA**

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

#### Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668 - 5324

## ©2012

Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. Form No. 11-246 111282

