

Model 270



The Gomco 270 is a tabletop gastric drainage aspirator for hospitals, clinics and home care. It has two pre-set drainage levels - 90mm Hg and 120 mm Hg - with 0.25 Ipm and 0.3 Ipm respective flow rates.

Key Features

- Silent operation due to a Gomco exclusive • thermotic pump design
- Compact, easy to store design takes up about one-half square feet of space for bedside use
- Floor bottle basket for stability with either bottle or canister
- Disposable hydrophobic bacteria filter helps prevent fluid and aerosol contamination

Specifications

Vacuum Setting: Flow Rate:

Sound Level: Pump Type: Electrical Cord: Dimensions:

Weight: UL/ETL Listed: Warranty: Regulatory Approvals: -90 mm Hg or -120 mm Hg 0.25 lpm @ -90 mm Hg and 0.3 @ -120 mm Hg Less than 47 dBA Thermotic 6' hospital grade 12" W x 8" D x 11" H (30 x 20 x 28 cm) 10 lbs. (5 kg) Yes 1 year FDA, HC

Model	Voltage	Glass Bottle	Disposable Canister
01-22-0270	115V/50-60 Hz	1 - 1200 ml	
01-22-2701	115V/50-60 Hz		1 - 1100 ml
01-73-0270	115V/50-60 Hz	1 - 1200 ml	

All models come with disposable suction tubing and 3 disposable hydrophobic bacteria filters.

01-90-2000	Disposable Suction Tubing Kit, one 18" & one 72" per pkg.	
01-90-3100	Disposable Hydrophobic Bacteria Filter, 1/8" NPT x 1/4" hose barb (3/pkg.)	
01-90-3489	Glass Collection Bottle, with cap & float assembly, 1200 ml	
01-90-3695	Disposable Collection Canister, stem inlet & hydrophobic filter, 1100 ml (12/case)	





Allied Healthcare Products, Inc. • 1720 Sublette Avenue • St. Louis, MO 63110 USA 314/771-2400 • 800/444-3954 • FAX: 314/771-0650 • www.alliedhpi.com

U.S. Toll Free: 800/444-3940

FAX: 800/477-7701

Canada & Mexico Toll Free: 800/446-0552 FAX: 800/246-6201

International Telephone: 314/268-1683 FAX: 314/771-5183

©Allied Healthcare Products, Inc.®

All specifications are nominal and subject to change without notice.

Form 21-00-0270 July 2013