

Lubricating Jelly











HR® Lubricating Jelly 10ml Pre-Filled Syringe is a sterile single use product, the oral tip design with easy-on & easy-off ergonomic cap is ideal for catheter use and medical kits.

Product Specs

MFR#	Item Description	Volume	HCPCS	Packaging
206-400	HR® Lubricating Jelly 10ml Pre-Filled Syringe - STERILE	10ml	A4332	400 ea/cs

Medical Applications

Colonoscopy

Enema Insertion

X-Ray

Urethral/Catheter

Pelvic Examination

Lithotripsy

Surgical Procedures

General Medical Use

Sigmoidoscopy

Ultrasound Procedures

Endotracheal tube insertion

H·R[®] Medical Supply

National Distribution Independent Distribution International Distribution



Features & Benefits

- Kosher Certified
- Water soluble
- Sterile Products
- Premium Viscosity
- Made in the USA
- Bacteriostatic
- Non-Spermicidal
- Chlorhexidine (CHG) Free
- Latex Free
- LUOM Labeling



Flip-Top Tubes





4oz. (113g) Flip-Top Tube

The 4 oz. Flip-Top Tube is our largest offering of sterile Lubricating Jelly. This size accommodates a broad range of medical procedures.

Product Specs

MFR#	Item Description	Volume	HCPCS	Packaging
201	HR® Lubricating Jelly 4oz. (113g) Flip-Top Tube - STERILE	4oz. (113g)	A4402	12 ea/bx - 6 bx/cs





20z. (56.5g) Flip-Top Tube

HR® Lubricating Jelly 2 oz. Flip-Top Tube is sterile and accommodates a broad range of medical procedures where a multi-use product is appropriate.

Product Specs

MFR #	Item Description	Volume	HCPCS	Packaging
203	HR® Lubricating Jelly 2oz. (56.5g) Flip-Top Tube - STERILE	2oz. (56.5g)	A4402	12 ea/bx - 8 bx/cs





20z. (56.5g) 2X Flip-Top Tube

HR® Lubricating Jelly 2 oz. 2X Flip-Top Tube is a sterile, multi-use option with a thicker consistency than the original Lubricating Jelly formula. The 2X formula is preferred by endoscopy physicians.

Product Specs

MFR#	Item Description	Volume	HCPCS	Packaging
210-192	HR® Lubricating Jelly 2X 2oz. (56.5g) Flip-Top Tube - STERILE	2oz. (56.5g)	A4402	192 ea/cs

Mission Statement

The mission of HR® is to be a supplier of personal care, consumer and medical supply products, consistently delivering superior product quality, extraordinary product performance and innovative ideas and products. We pledge to provide customized customer service and build the unconditional trust of our customers, so our customers can "Trust in Us®". Our mission is achievable by building a strong corporate culture, developed by our employees' commitments to our corporate values.







HR® Lubricating Jelly 3g OneShot® is a sterile, single use product designed to improve application, minimize waste, and significantly reduce the risk of cross contamination.

MFR#	Item Description	Volume	HCPCS	Packaging
207	HR® Lubricating Jelly 3g OneShot® - STERILE	3g	A4332	144 ea/bx - 6 bx/cs





5g OneShot®

HR® Lubricating Jelly 5g OneShot® is a sterile, single use product designed to improve application, minimize waste, and significantly reduce the risk of cross contamination.

MFR #	Item Description	Volume	HCPCS	Packaging
209	HR® Lubricating Jelly 5g OneShot® - STERILE	5g	A4332	144 ea/bx - 6 bx/cs





3g OneShot® CarePac™

The 3g OneShot® CarePac™ box of 30 is a convenient option for dispensing our sterile Lubricating Jelly to patients out of a physician's office/pharmacy or DME.

MFR #	Item Description	Volume	HCPCS	Packaging
208	HR [®] Lubricating Jelly 3g OneShot® CarePac™ - STERILE	3g	A4332	30 ea/bx - 48 bx/cs

Vision Statement

The vision of HR® is to be recognized as a global manufacturer and marketer of personal, consumer and medical products with a strong portfolio of globally recognized brands & products. This vision is achieved through the commitment of our employees to our shared mission and values.







1.25oz. (36g) OneShot®

HR® Lubricating Jelly 1.25 oz. OneShot® is a sterile single use alternative to our 2 oz. Flip-Top Tube. Research has shown that a 2 oz. tube is too much volume for some procedures. The 1.25 oz. OneShot® reduces waste and the risk of cross contamination.

Product Specs

MFR #	Item Description	Volume	HCPCS	Packaging
205	HR® Lubricating Jelly 1.25oz. (36g) OneShot® - STERILE	1.25oz. (36g)	A4402	48 ea/bx - 6 bx/cs





1.25oz. (36g) <mark>2X</mark> OneShot®

HR® Lubricating Jelly 2X OneShot® is a sterile, single use option with a thicker consistency than the original Lubricating Jelly formula. Our OneShot® is designed to significantly reduce cross contamination. The 2X formula is preferred by endoscopy physicians.

Product Specs

MFR#	Item Description	Volume	HCPCS	Packaging
211-288	HR® Lubricating Jelly 2X 1.25oz. (36g) OneShot® - STERILE	1.25oz. (36g)	A4402	288 ea/cs
211-576	HR® Lubricating Jelly 2X 1.25oz. (36g) OneShot® - STERILE	1.25oz. (36g)	A4402	576 ea/cs





1.25oz. (36g) OneShot® SafeWrap™

HR's SafeWrapTM packaging takes our OneShot® offering and makes it operating room and sterile-field friendly. This product is our original HR® Lubricating Jelly 1.25 oz. OneShot® packaged in a STERILE SafeWrapTM sleeve.

Product Specs

MFR#	Item Description	Volume	HCPCS	Packaging
205ST	HR [®] Lubricating Jelly 1.25oz. (36g) OneShot® SafeWrap™ - STERILE	1.25oz. (36g)	A4402	50 ea/bx - 4 bx/cs



HR Pharmaceuticals, Inc. 2600 Eastern Blvd., Suite 201 York, PA 17402





History


The legacy of HR® Lubricating Jelly dates back to 1925 when the HR® formula was imported from Germany. Soon thereafter, Holland-Rantos began marketing the original "Vaginal Jelly" to be known as HR® Lubricating Jelly. In the 1940's the Holland-Rantos Co., and its HR® branded product line was sold to Young's Rubber Co., which manufactured the Trojan® Condom.

With its rich history in hand, during the 1970's and 1980's, HR° Lubricating Jelly morphed itself into the go to medical lubricant that controlled a significant market share in the medical world. HR° Lubricating Jelly has always been considered a premium medical lubricant that offers UniversalTM usage to Physicians.

ENVIRONMENT POLICY

As a manufacturer of medical products we recognize that we have a responsibility to reduce our impact on the environment and to conserve natural resources. We endeavour to use the smallest viable packaging size to minimize material usage and reduce our landfill footprint.

When transporting products we make sure our operations are as efficient and environmentally friendly as possible.





Corporate Address: HR Pharmaceuticals, Inc.

2600 Eastern Blvd., Suite 201 York, PA 17402

www.HRpharma.com Office: 877-302-1110 Fax: 717-685-2590