

# Half Round, Stationary, Padded Phlebotomy Stand

#### **MODEL 6940**

- Vinyl upholstered, padded surface
- Easy-clean, 10-3/4" x 8-1/4" W x 3-1/2" D plastic bin (27cm L x 21cm W x 9cm D)
- Chrome plated steel pole with knob-lock adjustment
- 5-leg, cast aluminum base with black, non-marking, nylon glides

NOTE: Supplies pictured not included.

All blood drawing products are non-returnable unless they are unused in unopened, original packaging.



### **MODEL 6940 SPECIFICATIONS**

Тор	22" L x 11" Radius
Height Range	30" — 38"
Base	23" D



## **Prop 65 Warning**

**! WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

## **Upholstery Colors**



















Warm Gray