

The Handheld Pulse Oximeter is intended for continuous monitoring, spot-checking of functional pulse oxygen saturation (%SpO2) and pulse rate (PR) of adult and pediatric patients in nursing homes, hospitals and home care.

- Rubber grip design offers special protection
- · Compact, lightweight design for simple, one-hand operation
- High brightness LED displays SpO2, pulse rate and pulse bar
- Continuous operation for up to 30 hours
- · Meets the standard for conscious sedation
- Up to 99 patients' ID and 72-hour record storage
- Visual & audio alarm, low battery alarm
- Data transfer to PC for storage or printing
- Convenient 2 AA batteries (not included)
- Comes with one year limited warranty
- Available for purchase as 1 unit per case
- Not made with natural rubber latex



See our website for our Fingertip Pulse Oximeters

Finger Sensor

with cable

*Depending on market demand and Dynarex stock, we will be happy to address any orders for items currently not shown.













dynarex° is your one-stop medical supply company since 1967

For more information contact your local Dynarex® distributor sales representative. Or call **888-396-2739** to find your nearest distributor. 10 Glenshaw Street • Orangeburg, NY 10962 • www.dynarex.com

Packaged