

Hip Positioner & Pad

HCPCS: E1399-1,E1399



- Height adjustable seat harness positions hip and pelvis and allows user to shift weight from one leg to the other
- Shorten length of straps to assist in weight bearing or lengthen straps to use as a safety sling
- Position of straps on handlebar will affect users pelvis position
- Removable pad provides cushioned comfort

PRODUCTS

Description

TK 1070 L	Hip Positioner & Pad for Moxie, Lg,1/bx
TK 1070 S	Hip Positioner & Pad, Moxie, Sm, 1/bx