







Introduction

Chattanooga is the world's largest manufacturer of rehabilitation equipment for treating musculoskeletal, neurological and soft tissue disorders. For over six decades, Chattanooga has set a benchmark for leadership, reliability and excellence. Our products contribute to better treatment outcomes in hospitals, clinics and home settings worldwide. We lead the physiotherapy industry by example and, through continuous innovation, provide real world solutions for clinicians and their patients.

As the world's leading brand in rehabilitation, Chattanooga is sold in over 80 countries worldwide. Chattanooga is directly represented in our key markets but also has long lasting relationships with dedicated distributors in the USA, Canada, Latin America, EMEA and Asia Pacific. Our Chattanooga distributors and customers are supported for order placement, training, product delivery and after sales service by a dedicated team of highly skilled DJO employees based all over the world.

Chattanooga is a brand of DJO Global Inc. DJO Global provides solutions for musculoskeletal health, vascular health and pain management. Our products help prevent injuries or rehabilitate after surgery, injury or from degenerative disease, enabling patients to regain or maintain their natural motion.

Visit www.DJOglobal.eu





Contents

Introduction	1
Hydrocollator Mobile Heating Units	3
Hydrocollator Stationary Heating Units	4
Hydrocollator HotPacs	5
Terry Covers	6
Stainless Steel Carts	6
ColPac Chilling Unit	7
ColPac, boo-boo pac, Sport-Pac, Nylatex Wraps	8



Hot and Cold

Chattanooga is renowned for its "hot and cold products" all over the world

Chattanooga is the proud inventor of the original Hydrocollator HotPacs and heating unit, an all natural product that was first to market more than 60 years ago and is still the first choice of physiotherapists and other medical professionals when using heat therapy.

This HotPac generates deep moist heat which delivers increased advantages over dry heat to relieve chronic arthritic pain, acute muscle strains and general musculoskeletal discomfort.

As with moist heat therapy, Chattanooga sets the standard in cold therapy with the original ColPac range, which is now sold all over the world

3 Hot and Cold

Ordering Information		
Part Number	Description	
M-2 Mobile Hydro	ocollator	
2402	120 V Unit	
2407	220 V Unit	
Standard Accesso	ries	
1006	Standard size HotPacs (12x)	
M-4 Mobile Hydr	ocollator	
2502	120 V Unit	
2506	220 V Unit	
Standard Accesso	ries	
1006	Standard size HotPacs (24x)	
SS-2 Mobile Hydr	ocollator	
2302	120 V Unit	
2307	220 V Unit	
Standard Accessories		
1006	Standard size HotPacs (8x)	

Technical Specifications

M-2 Mobile Hydrocollator		
Mains Power:	110~240 V 50/60 Hz	
Power Consumption:	1000 W	
Weight:	39 kg	
Dimensions*:	67 x 40 x 84 cm	
Electrical Safety Class:	Class 1, Type B	
Safety Tests:	EN 60601-1	

M-4 Mobile Hydrocollator

Mains Power:	110~240 V 50/60 Hz
Power Consumption:	1500 W
Weight:	60 kg
Dimensions*:	89 x 51 x 84 cm
Electrical Safety Class:	Class 1, Type B
Safety Tests:	EN 60601-1

SS-2 Mobile Hydrocollator

Mains Power:	110~240 V 50/60 Hz
Power Consumption:	1000 W
Weight:	30 kg
Dimensions*:	53 x 41 x 84 cm
Electrical Safety Class:	Class 1, Type B
Safety Tests:	EN 60601-1

* Dimensions are listed as L x W x H

Hydrocollator[®] Mobile Heating Units





Features

Hydrocollator Mobile Heating Units

These durable stainless steel water tanks house the Hydrocollator Moist Heat HotPacs. The units are thermostatically controlled to ensure the ideal therapeutic temperature for HotPacs. Hydrocollator remains the benchmark heating unit against which all other units are judged.

- Full fiberglass insulation to prevent heat loss
- Dependable, rugged stainless steel design
- Simple to fill with water. No plumbing required
- Constant temperature of HotPacs is maintained
- Mobile units are equipped with 8 cm swivel, rubber casters for friction free movement about the clinic
- One year warranty

	M-2	M-4	SS-2	E-2	E-1
Tank Capacity	69 L	136 L	49 L	43 L	15 L
Temperature Range	71°-74°C	71°-74°C	71°-74°C	71°-74°C	71°-74°C
Thermal Cut-Out Temp.	82°-85°C	82°-85°C	82°-85°C	82°-85°C	82°-85°C
Accuracy	10%	10%	10%	10%	10%
Heat up Time (to 70°C)	8 Hours	8 Hours	4 Hours	4 Hours	2 Hours
Cool Down Time (from 70°C)	3 Hours	4 Hours	4 Hours	2 Hours	1 Hour
Fiberglass Insulation	Yes	Yes	Yes	No	No

www.chattgroup.eu 4

Hydrocollator[®] Stationary Heating Units

F. E-1

Hydrocollator Accessories

Side Table Rack

Stainless steel work surface and wire shelves add efficiency to units. Racks available to fit most models (M-4, M-2).

Lotion Warmers

Keeps massage lotion and ultrasound gel at a comfortable temperature. Easily attaches to any treatment table or wall. Mini-Lotion Warmer hangs on the front of the Hydrocollator E-1 or E-2 heating units (bottles not included).

Wall Mounted Towel Rack

This stainless steel rack provides a convenient way to dry or store towels and terry covers (towels not included).







Chattanooga®

Ordering Infor		
i ure itumber	Description	
E-1 Stationary Hydrocollator		
2102	120 V Unit	
2107	220 V Unit	
Standard Accesso	pries	
1006	Standard size HotPacs (4x)	
E-2 Stationary H	ydrocollator	
2802	120 V Unit	
2806	220 V Unit	
Standard Access	ories	
1006	Standard size HotPacs (3x)	
1004	Oversize HotPacs (2x)	
1002	Neck Contour HotPac (1x)	
Hydrocollator	Accessories	
Side Table Rack		
4010	Side Table Rack (not shown)	
4230	Extra Shelf for Side Table Rack	
TM-1 Electric Lot	ion Warmer	
4012	TM-1 Electric Lotion Warmer (110 V only)	
Mini Non-Electri	c Lotion Warmer	
22100	Mini Lotion Warmer	
Dial Thermomet	er	
4228	Dial Thermometer (not shown)	
Wall Mounted Towel Rack		
4016	Wall Mounted Towel Rack 5 x 84 cm	
Technical Spec	ifications	
E-1 Stationary Hy	/drocollator	
Mains Power:	110~240 V 50/60 Hz	
Power Consump	tion: 1000 W	
Weight:	7 kg	

Safety Tests: * Dimensions are listed as L x W x H

Dimensions*:

Safety Tests:

Mains Power:

Dimensions*:

Weight:

Power Consumption:

Electricial Safety Class:

Electricial Safety Class:

E-2 Stationary Hydrocollator

33 x 20 x 41 cm

Class 1, Type B

EN 60601-1

1000 W

38 x 33 x 51 cm

Class 1, Type B

EN 60601-1

14 kg

110~240 V 50/60 Hz

5 Hot and Cold

Ordering Information		
Part Number	Description	
Hydrocollator Ho	tPac - Moist Heat Hot Pack	
A 1008	23 cm x 39 cm	
B 1010	25 cm x 61 cm	
C 1004	Oversize 38 cm x 61 cm	
D 1002	Neck Contour 61 cm	
E 1014	Half Size 13 cm x 30 cm	
F 1006	Standard 25 cm x 30 cm	
G 1012	Knee-Shoulder 25 cm x 50 cm	

Hydrocollator[®] HotPacs[™]

Chattanooga®



Features

The Original Hydrocollator HotPac

Hydrocollator Moist Heat HotPacs are the first choice for physiotherapists, physicians, chiropractors, athletic trainers and kinesiologists for over 60 years.

Chattanooga is the inventor of the original Hydrocollator HotPac and proud to still be producing this all natural product today. The benefits of moist heat therapy over dry heat are well established for pain relief due to arthritis, myofibrosis, back pain, muscle strains and spasms.

Hydrocollator HotPacs

- Seven HotPac sizes designed to fit every body contour
- Provides up to 30 minutes of deep, soothing, moist heat
- Tabs on most HotPacs for easy rotation in busy clinics
- Cost effective, provides hundreds of treatments from one HotPac
- All natural pack containing bentonite, a natural hydrocolloid filling, wrapped in a woven cotton blend
- One year warranty



www.chattgroup.eu 6

Terry Covers



Features

All-Terry Covers

- For high-volume institutional use
- Hook and loop closures to prevent unwrapping or slipping

Foam-Filled Terry Covers

- For home or clinic use
- Polyurethane foam layer between two layers of terry
- Available with or without pockets

Stainless Steel Carts









chatt	anooga®	
Onon	anooga	

Ordering Information		
Part Number	Description	
All-Terry Covers		
1118	Standard 48 cm x 69 cm	
1126	Oversize 61 cm x 91 cm	
1104	Neck Contour 64 cm x 46 cm	
Foam-Filled Terry	Covers	
1102	Standard Size 48 cm x 69 cm	
1124	Oversize 61 cm x 91 cm	
1120	Neck Contour 64 cm x 41 cm	
Foam-Filled Pocke	t Terry Covers	
1108	Standard 43 cm x 69 cm	
1110	Oversize Pocket 61 cm x 91 cm	

Hydrocollator HotPac Sets (Combines HotPac and Foam Filled Terry Cover with retail packaging)

1064	Standard Set
1062	Neck Set

Stainless Steel Carts

4002	Model MB, Specially designed to hold the E-2 and E-1 Hydrocollator heating units 51 x 48 x 43 cm (20"H x 19"W x 17"D)
4004	Model MB-T, Accommodates stationary heating, ultrasound and electrotherapy units 76 x 48 x 43 cm (30"H x 19"W x 17"D)
4018	Model MB-TD, Accommodates stationary heating and electrotherapy units. Supply drawer included 76 x 48 x 43 cm (30"H x 19"W x 17"D)

7 Hot and Cold

Ordering Information		
Part Number	Description	
ColPac Chilling U	nit Model C-5	
3102	າາo V Unit	
3114	220 V Unit	
Standard Accessories		
1500	Standard Size ColPacs (6x)	
1506	Half Size ColPacs (6x)	

Technical Specifications

ColPac Chilling Unit	
Mains Power:	220~240 V 50 Hz; 115 V, 60 Hz
Weight:	74 kg
Dimensions*:	69 x 41 x 84 cm
Electricial Safety Class:	Туре В
Safety Tests:	EN 60601-1
* Dimensions are listed as L	хWхH

ColPac[®] Chilling Unit





Features

ColPac Chilling Unit

Large capacity and faster chilling time mean your patients will never have to wait for a ColPac. The ColPac Chilling Unit helps you keep pace throughout a busy day. No plumbing required, simply plug into a standard electrical outlet.

All ColPac Chilling Units include:

- Closed-cell foam insulation for efficiency
- Stainless steel, welded construction for long life
- Heavy-duty compressor for more efficient cooling
- Specially designed interior chilling coils cool packs faster
- 8 cm (3") swivel-type rubber casters for silent, friction-free movement of mobile units
- Drain valve for easy cleaning and defrosting
- One year warranty

www.chattgroup.eu 8

Colpac®, boo-boo pac™, Sport-Pac™, Nylatex® Wraps



Features

The Original ColPac

Cold therapy, when applied directly to the affected area, is used to control pain due to inflammation, swelling, multiple sclerosis, muscle strains, sprains and spasms.

- Seven different sizes of blue vinyl ensure effective cold therapy regardless of body area or shape. Also available are four sizes of black polyurethane
- Effective relief of pain
- Filled with a non-toxic silica gel that remains pliable throughout treatment
- Delivers up to 30 minutes of soothing relief from acute pain, swelling and fever
- Latex free
- One year warranty

boo-boo pac

- Fun, bear-shaped design especially made for kids
- Ideal for pediatric applications
- Filled with non-toxic silica gel
- Delivers up to 30 minutes of soothing relief from acute pain, swelling and fever

Nylatex Wraps

- Ideal for holding ColPacs in place
- Several sizes available
- Good for hundreds of uses

FlexiPac

- Microwaveable to provide soothing heat, or after freezer storage to provide penetrating cold therapy
- Flexible gel conforms to any body contour at any temperature

Sport-Pac

- Soccer ball design is ideal for youth sport applications
- Easily transferred from the freezer to the team cooler for sideline use
- Provides 30 minutes of effective cold therapy for acute pain due to strains, sprains and bruises

Chattanooga®

Ordering Information Part Number Description Blue Vinyl ColPac **A** 1512 Oversize 28 cm x 53 cm **B** 1500 Standard Size 28 cm x 36 cm **C** 1506 Half Size 19 cm x 28 cm **D** 1508 Neck contour 58 cm **E** 1502 8 cm x 28 cm **F** 1504 Quarter Size 14 cm x 19 cm

G 1510	Eye (Specially shaped to cover closed eyelids)
Black Polyurethane	e ColPac
H 1554	Neck Contour 53 cm Long
I 1556	Oversize 32 cm x 47 cm
J 1552	Standard Size 25 cm x 34 cm
K 1562	Half Size 17 cm x 28 cm

${\sf L}\,$ boo-boo pac (Shipped in a display box containing 12 packs)

1534 Royal Blue boo-b	000 pac
-----------------------	---------

M Nylatex Assortment (1 pack of 9)

1200	Nylatex Assortment (One of each size listed below except 6 cm x 61 cm
Nylatex Individua	I (1 pack of 3)
1202	6 x 46 cm
1204	6 x 91 cm
1206	6 x 122 cm
1208	10 x 46 cm
1210	10 X 91 CM
1212	10 X 122 CM
1213	10 X 152 CM
1214	15 x 46 cm
1216	15 x 91 cm
1218	15 X 122 CM
1219	15 X 152 CM
1264	6 х 61 cm

N FlexiPac Hot & Cold Compresses

4026	FlexiPac 3 cm x 15 cm (48/case)
4020	FlexiPac 13 cm x 25 cm (24/case)
4029	FlexiPac 20 cm x 35 cm (12/case)

• Sport-Pac (Shipped in a display box containing 10 packs)

Sport-Pac

1517

9 Hot and Cold	

DJO GLOBAL

AUSTRALIA: T: +1300 66 77 30 F: +1300 66 77 40 E: customerservice.au@DJOglobal.com

CHINA:

T: (8621) 6031 9989 F: (8621) 6031 9709 E: information_china@DJOglobal.com

GERMANY:

T: +49 761 4566 01 F: +49 761 456655 01 E: infoservice@DJOglobal.com

SOUTH AFRICA:

T: +27 (0) 87 3102480 F: +27 (0) 86 6098891 E: chattanooga.saf@DJOglobal.com

UK & IRELAND: T: +44 (0)1483 459 659 F: +44 (0)1483 459 470 E: ukorders@DJOglobal.com BENELUX: T: Belgium 0800 18 246 T: Netherlands 0800 0229442 T: Luxemburg 8002 27 42 E: benelux.orders@DJOglobal.com

DENMARK, FINLAND, NORWAY & SWEDEN: T: Sweden 040 39 40 00 T: Norway 8006 1052 T: Finland 0800 114 582 T: Denmark +46 40 39 40 00 E: info.nordic@DIOglobal.com

ITALY: T: +39 02 484 63386 F: +39 02 484 09217 E: vendite@DJOglobal.com

SPAIN: T: +34 934 803 202 **F:** +34 934 733 667 **E:** ventas@DJOglobal.com

UNITED STATES: T: +1 800 336 6569 F: +1 800 936 6569

E: customercare@DJOglobal.com

DJO GLOBAL, EXPORT CENTRES

ASIA-PACIFIC:

DJO Asia-Pacific Limited Unit 1905, 19/F, Tower II Grand Central Plaza 138 Shatin Rural Committee Road Shatin HONG KONG **T:** +852 3105 2237 **F:** +852 3105 1444 **E:** info.asia@DJOglobal.com EUROPE, MIDDLE EAST & AFRICA: DJO Benelux Welvaartstraat 8 2200 Herentals BELGIUM

E: info.emea@DJOglobal.com

CANADA: T: +1 1866 866 5031 **F:** +1 1866 866 5032 **E:** canada.orders@DJOglobal.com

FRANCE: T: +33 (0)5 59 52 86 90 F: +33 (0)5 59 52 86 91 E: sce.cial@DJOglobal.com

INDIA: T: +91 44 66915127 E: customercare.india@DJOglobal.com

SWITZERLAND: T: +41 (0) 21 695 2360 F: +41 (0) 21 695 2361 E: info@compex.ch

LATIN AMERICA:

DJO Global, Inc 1430 Decision Street Vista CA 92081-8553 U.S.A. **T:** 1 800 336 6569 **F:** 1 800 936 6569 **E:** info.latam@DJOglobal.com



DJO Global | 1a Guildford Business Park | Guildford | Surrey | GU2 8XG | UNITED KINGDOM www.DJOglobal.eu