

## Hydraulic, Deluxe Silver Vein Patient Lift

HCPCS: E0630



- Silver Vein steel construction
- 6 point Swivel Bar that can accommodate 2 or 4 sling straps plus 2 chain connections
- High-performance hydraulics raise or lower individuals gradually and safely from any stationary position
- Choose either 3" or 5" casters
- Adjustable-width base

PRODUCTS	Description
13023SV	Patient Lift,New Style,Silver Vein,1/cs

SPECIFICATIONS	Description
Base Height	6.5" (13023SV)
Base Width	42" (open) / 22.5" (closed)
Boom Height	17.5" - 72" (13023SV)
Carton Shipping Weight	84 lbs.
Casters	5" (13023SV)
Warranty	Limited Lifetime (frame), 1 Year Limited (pump)
Weight	72 lbs.
Weight Capacity	450 lbs

PARTS	Description
13023SV01	HangerAttachmentBolt13023SV,1ea, EA
13023SV03	MastComplete13023SV,1ea, EA
13023SV04	UpperMast13023SV,1ea, EA
13023SV06	PumpAttachmentBolt13023SV,1ea, EA
13023SV07	MastHingeBolt13023SV,1ea, EA
13023SV08	PushHandleGrip13023SV,1PR, PR
13023SV09	Pump13023SV,1ea, EA
13023SV09-FPHH	Pump Handle Hardware, 13023SV, 1/ea, EA
13023SV10	SpreaderBarHandleGrip13023SV,1ea, EA
13023SV11	SpreaderBar13023SV,1ea, EA
13023SV12	ThumbScrew13023SV,1ea, EA
13023SV13	SpreaderBarReceiver13023SV,1ea, EA
13023SV16	RearCasterWithLock13023SV,1ea, EA
13023SV18	FrontCasterW/oLock13023SV,1ea, EA
13023SV19	BaseTipCover13023SV,1ea, EA

**PARTS****Description**

13023SV20	TubePlug13023SV,1ea, EA
13023SV21	CradlePivot13023SV,1ea, EA
13023SV22	PumpHandGrip13023SV,1ea, EA
13023SV23	MastTipCover13023SV,1ea, EA
13052RL	ReleaseLever,13052,1/ea, EA