

hand dynamometers

Baseline[®] digital smedley dynamometer

- digital, lightweight smedley style dynamometer offers accurate grip strength results
- 5-preset grip sizes
- captures and displays maximum grip force
- 200 lb (90 kg) capacity
- assesses results according to age and gender
- saves/stores results for up to 19 users for fast retrieval
- requires 2 AAA batteries included



12-0286 digital smedley 110.00

Baseline[®] smedley spring dynamometer

- adjustable handle comfortably accommodates any hand size
- dial indicator remains at maximum until reset
- calibrated in both lb and kg
- includes protective carrying case



12-0281 adult (220 lb/100kg) 200.00
12-0282 child (110 lb/50kg) 200.00

Baseline[®] pneumatic squeeze dynamometers

single bulb three-bulb set



- the inexpensive way to measure grip and pinch strength
- available with and without a maximum force indicator (reset) that remains at the maximum reading until reset
- 2" bulb, fluid-filled
- calibrated in PSI
- measures both grip and pinch strength
- similar to the Martin Vigorometer
- three different sized bulbs (1", 1½" and 2") can be fitted to the gauge
- lightweight and portable
- maximum reading remains until reset
- requires no fluid
- unit is calibrated in PSI (0-30) and kPa (0-200)

with reset
12-0291 30 psi 75.00
12-0293 15 psi 75.00

without reset
12-0290 30 psi 50.00 12-0297 3-bulb dynamometer 150.00
12-0292 15 psi 50.00

Jamar[®] dynamometers

hydraulic digital Jamar+[®]



- hydraulic dynamometer gives accurate and repeatable grip strength readings
- adjustable 5-position handle
- maximum strength indicator remains until reset
- reads lb and kg
- 200 lb (90 kg) capacity
- for routine screening or grip strength evaluation
- sturdy aluminum body construction
- rapid exchange test with audible signal
- easy-to-read LCD display can be set to display lb or kg
- 200 lb (90 kg) capacity
- allows up to 5 trials for left and right hand
- auto calculates average standard variation and coefficient of variation

12-0600 Jamar[®] hydraulic 365.00 12-0604 Jamar+[®] digital 365.00

evaluation