


## Economy 5-Leg, 4-Hook IV Pole



### MODEL: IV-474

- 4-hook top
  - Chrome plated steel pole/twist-lock adjustment
  - 5-leg, nylon composite base
  - Detachable ram's horn hooks
  - Dual wheel casters
  - See available options on next page
- 45 lbs. load capacity under normal use


*Shipping Weight: 15 lbs.*

#### SPECIFICATIONS

Base Diameter—21"

Height Range: 46"—83"

**California residents please be advised of the following, pursuant to Proposition 65:**

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1

**OPTIONS:**



*IV-46 Pump Tray Support*



*IV-41 Universal Clamp*



*IV-55 Extra Hooks and Organizer Hooks*



*Heavy Duty Black Epoxy Wire Basket*  
*IV-51B - 14" x 8" x 8.5"*  
*IV-52B - 12" x 6" x 6.5"*



*Heavy Duty Stainless Steel Wire Basket*  
*IV-51S - 14" x 8" x 8.5"*  
*IV-52S - 12" x 6" x 6.5"*



*IV-57 Single Ram's Horn*



*IV-44 Oxygen Cylinder Holder*



*IV-43 Steering Handle*



*IV-56 Drainage Hook*