

## Stainless Steel Replacement Trays



### MODEL: M-24

- Rimmed, stainless steel tray

\*Clinton stainless steel items are not M.R.I. compatible.


*Shipping Weight: 5 lbs.*

#### SPECIFICATIONS

Base: 18" W x 22" D

Height Range: 31"—50"

California residents please be advised of the following, pursuant to Proposition 65:

 **WARNING:** This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

1