

Integra[®] Miltex[®]
Dermatology Instruments



Integra facility
Rietheim-Weilheim, Germany





Limit uncertainty with quality instruments from the brand you trust

For more than 50 years, Integra Miltex has established itself as one of the largest and most respected suppliers of hand-held surgical instruments. Today, Integra Miltex operates as a business unit of Integra LifeSciences Holdings Corporation (“Integra”). Integra is a diversified medical technology company that develops, manufactures, and markets medical devices for use in a variety of applications.

Integra Miltex, through Miltex GmbH in Rietheim-Weilheim, Germany, continues to employ the skills of the finest instrument craftsmen in the world. The surgical ‘feel’ of our instruments, superb cutting ability, and smoothly beveled box locks are just a few features that define our quality and performance.


This brochure represents a comprehensive offering of specialty instruments from Integra for the dermatologist. The selection includes premium SuperCut, Tungsten Carbide and Ceramic Coated scissors as well as a variety of forceps, needle holders, dermal curettes, hooks, extractors, disposable products, and care & handling solutions.

Integra® Miltex®

Specialty Scissors


Integra Miltex SuperCut Scissors

SuperCut scissors feature a razor-sharp upper blade and a lower blade with micro-fine serrations to hold tissue in place and help prevent slippage. Integra Miltex SuperCut scissors have two black ring handles for easy identification.



Integra Miltex Tungsten Carbide Scissors

These premium, stainless steel scissors, identified by two gold ring handles, have tungsten carbide inserts seamlessly bonded to the blades. Tungsten carbide scissors offer precision cutting performance and durability.





Integra Miltex SuperCut/Tungsten Carbide Scissors

Durability and precision cutting performance come together without compromise in SuperCut/Tungsten Carbide scissors, identified by one short and one long gold ring handle.

Tungsten Carbide insert honed to a razor-sharp edge

Tungsten Carbide insert with micro-fine serrations




Integra® Padgett® Ceramic Coated Scissors

Padgett ceramic scissors offer a threefold design that combines our SuperCut blade design with ceramic and tungsten carbide materials. The result is lasting cutting performance - approximately 81% longer than standard stainless steel scissors.

Tungsten Carbide insert honed to a razor-sharp edge

Tungsten Carbide insert with micro-fine serrations




Scissors

Stainless Steel, SuperCut, Tungsten Carbide,
SuperCut/Tungsten Carbide & Ceramic Coated



5-300



18-1453



18-SC-1652



5-SC-304TC



PM-5SC306TC

Iris Scissors

Reference	Description	Length
5-300	Straight	4" (102 mm)
5-302	Curved	4" (102 mm)
5-304	Straight	4½" (114 mm)
5-SC-304	Straight, SuperCut	4½" (114 mm)
5-304TC	Straight, Tungsten Carbide	4½" (114 mm)
5-SC-304TC	Straight, SuperCut/Tungsten Carbide	4½" (114 mm)
PM-5SC304CTC	Straight, SuperCut/Tungsten Carbide/Ceramic Coated	4½" (114 mm)
5-306	Curved	4½" (114 mm)
5-SC-306	Curved, SuperCut	4½" (114 mm)
5-306TC	Curved, Tungsten Carbide	4½" (114 mm)
5-SC-306TC	Curved, SuperCut/Tungsten Carbide	4½" (114 mm)
PM-5SC306CTC	Curved, SuperCut/Tungsten Carbide/Ceramic Coated	4½" (114 mm)

Gradle Scissors

Reference	Description	Length
18-1652	Slightly Curved	3¾" (95 mm)
18-SC-1652	Slightly Curved, SuperCut	3¾" (95 mm)
18-1653	Slightly Curved, Ribbon-Type Ring Handles	4" (102 mm)

Ragnell Dissecting Scissors

Reference	Description	Length
5-290	Curved	5" (127 mm)
5-SC-290	Curved, SuperCut	5" (127 mm)

Strabismus Scissors

Reference	Description	Length
5-312	Straight	4" (102 mm)
5-314	Curved	4" (102 mm)
5-312TC	Straight, Tungsten Carbide	4½" (114 mm)
5-314TC	Curved, Tungsten Carbide	4½" (114 mm)
18-1453	Straight, Ribbon-Type Ring Handles	4¼" (108 mm)
18-1454	Curved, Ribbon-Type Ring Handles	4⅛" (105 mm)

Scissors

Stainless Steel, SuperCut,
Tungsten Carbide & Ceramic Coated





PM-5SC18oCTC

Stevens Tenotomy Scissors

Reference	Description	Length
18-1460	Straight, Short Blades, Sharp	4 1/8" (105 mm)
18-1464	Curved, Short Blades, Sharp	4 1/8" (105 mm)
18-SC-1472	Straight, Long Blades, Blunt, SuperCut	4 1/2" (114 mm)
18-SC-1476	Curved, Long Blades, Blunt, SuperCut	4 1/2" (114 mm)

Wilmer Converse Scissors

Reference	Description	Length
18-1440	Blades Angled on Flat	4" (102 mm)
18-1440TC	Blades Angled on Flat, Tungsten Carbide	4 1/2" (114 mm)

Metzenbaum Scissors

Reference	Description	Length
5-179	Straight	5 1/2" (140 mm)
5-SC-179	Straight, SuperCut	5 1/2" (140 mm)
5-180	Curved	5 1/2" (140 mm)
5-SC-180	Curved, SuperCut	5 1/2" (140 mm)
5-180TC	Curved, Tungsten Carbide	5 1/2" (140 mm)
PM-5SC18oCTC	Curved, SuperCut/Tungsten Carbide/Ceramic Coated	5 1/2" (140 mm)
5-181	Straight	7" (178 mm)
5-182	Curved	7" (178 mm)
5-SC-182	Curved, SuperCut	7" (178 mm)
5-182TC	Curved, Tungsten Carbide	7" (178 mm)
PM-5SC182CTC	Curved, SuperCut/Tungsten Carbide/Ceramic Coated	7" (178 mm)

O.R. Scissors

Reference	Description	Length
5-2	Straight, Sharp/Sharp	4 1/2" (114 mm)
5-12	Straight, Sharp/Blunt	4 1/2" (114 mm)

Stitch Scissors

Reference	Description	Length
9-100	Spencer, Straight, Delicate, with Hook	3 1/2" (89 mm)
9-101	Shortbent, Curved, Delicate, with Hook	3 1/2" (89 mm)
9-102	Littauer Jr., Straight, with Hook	4 1/2" (114 mm)

Forceps

Stainless Steel & Tungsten Carbide



6-120



18-868



18-786



18-1107



7-24



6-118TC

Adson Dressing Forceps

Reference	Description	Length
6-118	Serrated, Delicate	4 ³ / ₄ " (121 mm)
6-118XL	Serrated, Delicate, Fenestrated Handles	4 ³ / ₄ " (121 mm)
6-118TC	Cross-Serrated, Delicate, Tungsten Carbide	4 ³ / ₄ " (121 mm)
17-2510	Serrated, Extra Delicate, 0.6 mm	4 ³ / ₄ " (121 mm)

Adson Tissue Forceps

Reference	Description	Length
6-120	1 x 2 Teeth, Delicate	4 ³ / ₄ " (121 mm)
6-120XL	1 x 2 Teeth, Fenestrated Handles	4 ³ / ₄ " (121 mm)
6-121	1 x 2 Teeth, Cross-Serrated, Delicate	4 ³ / ₄ " (121 mm)
6-123	1 x 2 Teeth, Tying Platform	4 ³ / ₄ " (121 mm)
6-123XL	1 x 2 Teeth, Tying Platform, Fenestrated Handles	4 ³ / ₄ " (121 mm)
6-123TC	1 x 2 Teeth, Cross Serrated Tying Platform, Tungsten Carbide	4 ³ / ₄ " (121 mm)
17-2500	1 x 2 Teeth, Extra Delicate, 0.5 mm	4 ³ / ₄ " (121 mm)

Brown-Adson Forceps

Reference	Description	Length
6-124	7 x 7 Teeth	4 ³ / ₄ " (121 mm)
6-124XL	7 x 7 Teeth, Fenestrated Handles	4 ³ / ₄ " (121 mm)

Bishop-Harmon Forceps

Reference	Description	Length
18-867	1 x 2 Teeth, Delicate, 0.3 mm	3 ³ / ₈ " (86 mm)
PM-4792	1 x 2 Teeth, Delicate, 0.4 mm	3 ¹ / ₄ " (83 mm)
PM-4791	1 x 2 Teeth, Standard, 0.6 mm	3 ¹ / ₄ " (83 mm)
18-868	1 x 2 Teeth, Standard, 0.7 mm	3 ³ / ₈ " (86 mm)

Iris Tissue Forceps

Reference	Description	Length
18-786	1 x 2 Teeth, Straight, 0.8 mm	4" (102 mm)

Swiss Cilia Forceps


Reference	Description	Length
18-1107	Diagonal Jaws, 3 mm	3 ³ / ₄ " (95 mm)


Ring-Handled Forceps

Reference	Description	Length
7-24	Hartman Mosquito Forceps, Straight	4" (102 mm)
7-26	Hartman Mosquito Forceps, Curved	4" (102 mm)
7-2	Halsted Mosquito Forceps, Straight	5" (127 mm)
7-4	Halsted Mosquito Forceps, Curved	5" (127 mm)
7-504	Backhaus Towel Clamp	3 ¹ / ₂ " (89 mm)
7-506	Backhaus Towel Clamp	5 ¹ / ₄ " (133 mm)

Needle Holders

Stainless Steel & Tungsten Carbide





18-1828

Webster Needle Holders

Reference	Description	Length
8-6	Smooth Jaws, Delicate	5" (127 mm)
8-6TC	Smooth Jaws, Delicate, Tungsten Carbide	5 ¹ / ₄ " (133 mm)
8-7	Smooth Jaws, Extra Delicate	4 ³ / ₄ " (121 mm)
8-7TC	Smooth Jaws, Extra Delicate, Tungsten Carbide	4 ³ / ₄ " (121 mm)

Halsey Needle Holders

Reference	Description	Length
8-8	Smooth Jaws	5" (127 mm)
8-8TC	Smooth Jaws, Tungsten Carbide	5" (127 mm)
8-8A-TC	Serrated Jaws, 4,000 Teeth Per Square Inch, Tungsten Carbide	5" (127 mm)

Baumgartner Needle Holders

Reference	Description	Length
8-40	Serrated Jaws	5 ¹ / ₂ " (140 mm)
8-40TC	Serrated Jaws, 4,000 Teeth Per Square Inch, Tungsten Carbide	5 ¹ / ₂ " (140 mm)

Mayo-Hegar Needle Holders

Reference	Description	Length
8-42	Serrated Jaws	5" (127 mm)
8-42TC	Serrated Jaws, Tungsten Carbide	5" (127 mm)
8-44	Serrated Jaws	6" (152 mm)
8-44TC	Serrated Jaws, Tungsten Carbide	6 ¹ / ₄ " (159 mm)

Derf Needle Holders

Reference	Description	Length
8-4	Horizontal Serrations	4 ³ / ₄ " (121 mm)
8-4TC	Serrated Jaws, 4,000 Teeth Per Square Inch, Tungsten Carbide	4 ³ / ₄ " (121 mm)

Olsen-Hegar Needle Holders with Suture Scissors

Reference	Description	Length
8-14	Serrated Jaws, Extra Delicate	4 ³ / ₄ " (121 mm)
8-14TC	Serrated Jaws, Extra Delicate, Tungsten Carbide	4 ³ / ₄ " (121 mm)
8-15	Serrated Jaws	5 ¹ / ₂ " (140 mm)
8-15TC	Serrated Jaws, Tungsten Carbide	5 ¹ / ₂ " (140 mm)

Castroviejo Needle Holders

Reference	Description	Length
18-1828	Smooth Jaws, Straight, with Lock	5 ¹ / ₂ " (140 mm)
18-1832	Smooth Jaws, Curved, with Lock	5 ¹ / ₂ " (140 mm)

Dermal Currettes

Stainless Steel





33-15

Fox Dermal Curettes

Reference	Description	Tip Diameter	Length
33-5B	Round	1 mm	5½" (140 mm)
33-5A	Round	2 mm	5½" (140 mm)
33-5	Round	3 mm	5½" (140 mm)
33-4	Round	4 mm	5½" (140 mm)
33-3	Round	5 mm	5½" (140 mm)
33-2	Round	6 mm	5½" (140 mm)

Piffard Dermal Curettes

Reference	Description	Tip Diameter	Length
33-6	Oval	Size 1 (3 mm)	6¼" (159 mm)
33-8	Oval	Size 2 (5 mm)	6¼" (159 mm)
33-10	Oval	Size 3 (6 mm)	6¼" (159 mm)
33-11	Oval	Size 4 (8 mm)	6¼" (159 mm)
33-11A	Oval	Size 5 (10 mm)	6¼" (159 mm)

Piffard Dermal Curettes, Narrow Handle

Reference	Description	Tip Diameter	Length
33-16-00	Oval	Size 00 (1 mm)	5½" (140 mm)
33-16-0	Oval	Size 0 (1.5 mm)	5½" (140 mm)
33-16-1	Oval	Size 1 (3 mm)	5½" (140 mm)
33-16-2	Oval	Size 2 (5 mm)	5½" (140 mm)
33-16-3	Oval	Size 3 (6 mm)	5½" (140 mm)

Miltex Dermal Curettes

Reference	Description	Tip Diameter	Length
33-12	Oval	Spoon Size 0 (3 mm)	6¼" (159 mm)
33-13	Oval	Spoon Size 1 (4 mm)	6¼" (159 mm)
33-14	Oval	Spoon Size 2 (5 mm)	6¼" (159 mm)
33-15	Oval	Spoon Size 3 (6 mm)	6¼" (159 mm)

Extractors & Dermal Punches

Stainless Steel



33-200

33-201

33-202


33-203

33-214

33-206

33-215

33-22



33-30

Schamberg Comedone Extractors

Reference	Description	Length
33-200	Square Loop Ends	3¾" (95 mm)
33-201	Standard Pattern with Crimped Small Loop	3¾" (95 mm)
33-202	Standard Pattern with Extra Fine Crimped Loop	3¾" (95 mm)

Miltex Comedone Extractors

Reference	Description	Length
33-203	Double-Ended, Flat and Wire End Expressors	4" (102 mm)
33-214	Double-Ended, One Fenestrated Cup and One Delicate, Curved Lancet with Sheath	6" (152 mm)

Saalfield Comedone Extractors

Reference	Description	Length
33-206	Extractor with Acne Lancet	4¼" (108 mm)
33-204	Extractor with Acne Lancet	5¾" (146 mm)

Unna Comedone Extractor

Reference	Description	Length
33-210	Spoon Ends	5¾" (146 mm)

Greg Comedone Extractor

Reference	Description	Length
33-215	One Fenestrated Cup and Chuck with Interchangeable No. 65 Lancet Blade	5" (127 mm)

Keyes Dermal Punches

Reference	Description	Length
33-20	2 mm	4" (102 mm)
33-22	3 mm	4" (102 mm)
33-24	4 mm	4" (102 mm)
33-25	5 mm	4" (102 mm)
33-26	6 mm	4" (102 mm)
33-28	8 mm	4" (102 mm)
33-30	Handle with Screw-In Keyes Knives: 2, 3, 4, 5, 6, 7 mm Mounted on a Storage and Sterilization Rack	3½" (89 mm)

Skin Hooks, Retractors & Nail Nippers

Stainless Steel





40-219

Joseph Skin Hooks

Reference	Description	Length
21-153	Single Prong, Sharp	6¼" (159 mm)
21-154	Double Prong, Sharp, 2 mm Wide	6¼" (159 mm)
21-156	Double Prong, Sharp, 5 mm Wide	6¼" (159 mm)
21-158	Double Prong, Sharp, 7 mm Wide	6¼" (159 mm)
21-160	Double Prong, Sharp, 10 mm Wide	6¼" (159 mm)

Miltex Skin Hook

Reference	Description	Length
21-152	Single Prong, Sharp, 2 mm	4¾" (121 mm)

Freer Skin Hooks

Reference	Description	Length
21-102	Double Prong, Sharp, 2.3 mm Wide	6" (152 mm)
21-103	Double Prong, Blunt, 3 mm Wide	6" (152 mm)

Frazier Skin Hooks

Reference	Description	Length
26-1084	Single Prong, Sharp, 5 mm	5" (127 mm)
26-1086	Single Prong, Blunt, 5 mm	5" (127 mm)

Miltex Cross Action Retractors

Reference	Description	Length
11-14	4 x 4 Sharp Prongs	4" (102 mm)
11-16	4 x 4 Blunt Prongs	4" (102 mm)

Desmarres Retractors

Reference	Description	Length
18-110	Size 1 (12 mm)	5½" (140 mm)
18-112	Size 2 (14 mm)	5½" (140 mm)

Nail Nippers and Nail Splitter

Reference	Description	Length
40-212-SS	Straight Jaws, Double Spring	5½" (140 mm)
40-219	Concave Jaws, Double Action	6" (152 mm)
40-230	English Anvil Nail Splitter, Double Spring	5" (127 mm)

Mohs Surgery


Instrument Set Up - A¹

Reference	Description	Length
33-5	Fox Dermal Curette, Round, 3mm Diameter	5½" (140 mm)
18-868	Bishop-Harmon Forceps, 1 x 2 Teeth, Standard, 0.7 mm	3¾" (86 mm)
4-7	Scalpel Handle No. 3, Stainless Steel, Graduated in mm and cm	5" (127 mm)
4-115	Surgical Blades #15, Carbon Steel, Sterile, Single Use, 100/Box	-
6-120	Adson Tissue Forceps, 1 x 2 Teeth, Delicate	4¾" (121 mm)
5-304	Iris Scissors, Straight	4½" (114 mm)
6-118	Adson Dressing Forceps, Serrated, Delicate	4¾" (121 mm)
4-100	Blade Safe Remover, Compatible with All Handle Styles	-



Instrument Set Up - A¹
Scan here for more information





Instrument Set Up - B¹


Reference	Description	Length
7-4	Halsted Mosquito Forceps, Curved	5" (127 mm)
6-123	Adson Tissue Forceps, 1 x 2 Teeth, Tying Platform	4¾" (121 mm)
5-2	O.R. Scissors, Straight, Sharp/Sharp	4½" (114 mm)
5-SC-290	Ragnell Dissecting Scissors, Curved, SuperCut	5" (127 mm)
8-42TC	Mayo-Hegar Needle Holder, Serrated Jaws, Tungsten Carbide	5" (127 mm)



Instrument Set Up - B¹
Scan here for more information

Disposable Biopsy Punches, Curettes & Biopsy Blade

Sterile, Single Use





33-100

Standard Biopsy Punches - 50/box

Reference	Size	Diameter (Sized to Scale)
33-31AA	1.0 mm	◦
33-31A	1.5 mm	◦
33-31	2.0 mm	◦
33-31B	2.5 mm	○
33-32	3.0 mm	○
33-33	3.5 mm	○
33-34	4.0 mm	○
33-35	5.0 mm	○
33-36	6.0 mm	○
33-37	8.0 mm	○
33-38	Assorted Sizes	

Biopsy Punches with Plunger - 25/box

Reference	Size	Diameter (Sized to Scale)
33-31AA-P/25	1.0 mm w/plunger	◦
33-31A-P/25	1.5 mm w/plunger	◦
33-31-P/25	2.0 mm w/plunger	◦
33-32-P/25	3.0 mm w/plunger	○
33-34-P/25	4.0 mm w/plunger	○

Dermal Curettes - 50/box


Reference	Size	Diameter (Sized to Scale)
33-52	2.0 mm	○
33-53	3.0 mm	○
33-54	4.0 mm	○
33-55	5.0 mm	○
33-57	7.0 mm	○

Miltex Biopblade® - 50/box

Reference	Description
33-100	Flexible Biopsy Blade

Disposable Scalpels

Sterile, Single Use





4-515C


4-522

Standard Scalpels - 10/box

Reference	Size
4-410	10
4-411	11
4-412	12
4-415	15
4-420	20
4-421	21
4-422	22
4-423	23

Safety Scalpels - 10/box

Reference	Size
4-510	10
4-511	11
4-512	12
4-515	15
4-515C	15C Angled Blade
4-522	22



** Sized to Scale

Scalpel Handles & Surgical Blades





** Sized to Scale

Reusable Scalpel Handles

Reference	Description	Size	Fits Blades 10-15C	Fits Blades 20-25
4-7	No. 3, Stainless Steel, Graduated in mm and cm	5" (127 mm)	✓	-
4-8	No. 4, Stainless Steel, Graduated in mm and cm	5 3/8" (136 mm)	-	✓
4-9	No. 7, Stainless Steel	6 1/2" (165 mm)	✓	-
4-10	No. 3L, Stainless Steel	8 1/2" (216 mm)	✓	-
4-16	No. 9, Stainless Steel	4 1/2" (114 mm)	✓	-
4-17	Siegel, Stainless Steel, Round Handle	5 3/4" (146 mm)	✓	-
4-400*	Miniature Blade Handle w/Chuck, Round Knurled (3K Type)	3 3/4" (95 mm)	-	-
4-401*	Miniature Blade Handle w/ Self-Locking Chuck, Round Knurled (3K Type)	3 3/4" (95 mm)	-	-
4-405*	Chisel Blade Handle, Self-Locking	5 1/4" (133 mm)	-	-

* Fits Mini Blade (blades not available from Integra Miltex)

Surgical Blades, Sterile - 100/box

Carbon Steel	Stainless Steel	Silicone Coated	Size
4-110	4-310	4610	10
4-111	4-311	4611	11
4-112	4-312	-	12
4-112B ^Δ	4-312B ^Δ	-	12B ^Δ
4-115	4-315	4615	15
4-115C ^{ΔΔ}	4-315C ^{ΔΔ}	-	15C ^{ΔΔ}
4-120	4-320	-	20
4-121	4-321	-	21
4-122	4-322	-	22
4-123	4-323	-	23
4-125	4-325	-	25

B^Δ is a double edged blade C^{ΔΔ} is an angled blade

Surgical Blade Remover

Reference	Description
4-100	Blade Safe Remover, Compatible with All Handle Styles

Instrument Care & Cleaning Products





3-800

Instrument Prep Enzyme Foam

Reference	Description
3-760	24 oz. (0.71 liter), Spray Bottle, Case of 12

EZ-Zyme® All Purpose Enzyme Cleaner

Reference	Description
3-750	3/4 oz. (22ml) Packets, Dispenser Box of 32, Contents of 1 Box makes 32 Gallons of Solution
3-755	1 Gallon (3.79 liter), Contents of 1 Gallon Container makes 170 Gallons of Solution

Instrument Cleaner

Reference	Description
3-720	8 fl. oz. (0.24 liter) Bottles, Case of 12, Contents of one 8 oz. Bottle makes 32 Gallons of Solution
3-725	1 Gallon Bottle (3.79 liter), Contents of 1 Gallon Container makes 512 Gallons of Solution

Instrument Cleaner & Lubricant

Reference	Description
3-710	1 Gallon (128 fl. oz.)

Spray Lube

Reference	Description
3-700	8 fl. oz. (0.24 liter) Bottles, Case of 12

Instrument Lubricant

Reference	Description
3-705	1 Gallon (128 fl. oz.)

Instrument Stain Remover

Reference	Description
3-740	3 oz. (85 grams), Case of 12

Nylon & Stainless Steel Brushes

Reference	Description
3-1000	Nylon Bristles, Set of 3, 7¼" (184 mm)
3-1001	Stainless Steel Bristles, Set of 3, 7¼" (184 mm)

Instrument Care Kit

Reference	Description
3-800	Instrument Cleaner (3-720), Spray Lube (3-700), Instrument Stain Remover (3-740), Instrument Care DVD (3-810), Nylon Brush and Stainless Steel Brush


Medicine Cup, Trays & Cassettes



3-914



3-921



4-6835



3-200200



3-939



3-938

Medicine Cup

Reference	Description	Dimensions
3-914	Graduated	2 ⁵ / ₃₂ " x 2 ⁵ / ₃₂ " (2 oz.)

Mayo Trays

Reference	Description	Dimensions
3-921	Perforated Tray	10" x 6 ¹ / ₂ " x 23/ ₃₂ " (254 mm x 165 mm x 18 mm)
3-922	Perforated Tray	14" x 10" x 9/ ₁₆ " (356 mm x 254 mm x 16 mm)
3-924	Perforated Tray	17" x 11 ⁵ / ₁₆ " x 3/ ₄ " (432 mm x 295 mm x 19 mm)
3-925	Perforated Tray	19 ¹ / ₁₆ " x 12 ²³ / ₃₂ " x 33/ ₆₄ " (487 mm x 323 mm x 13 mm)

Universal Cassette & Sterilization Trays

Reference	Description	Dimensions
4-6835	Universal Cassette without Mat	6 ¹ / ₂ " x 6 ¹ / ₄ " x 1 ³ / ₁₆ " (165 mm x 159 mm x 30 mm)
4-6835MAT	Silicone Mat with Holes, for Cassette 4-6835	6 ¹ / ₂ " x 6 ¹ / ₄ " (165 mm x 159 mm)
3-200100	Instrument Sterilization Tray with Silicone Mat, Single Layer	7" x 11" (178 mm x 279 mm)
3-200200	Instrument Sterilization Tray with Silicone Mat, Double Layer	7" x 11" (178 mm x 279 mm)
3-200300	Instrument Sterilization Tray with Silicone Mat, Single Layer	4" x 8" (102 mm x 203 mm)

Cath Trays & Covers

Reference	Description	Dimensions
3-937	Cath Tray, 8 x 5 x 1	8 ²⁵ / ₃₂ " x 4 ⁶¹ / ₆₄ " x 1 ⁴⁹ / ₆₄ " (223 mm x 126 mm x 45 mm)
3-938	Cath Tray, 8 ⁷ / ₈ x 5 x 3	8 ⁷ / ₈ " x 5" x 3" (225 mm x 127 mm x 76 mm)
3-939	Flat Cath Cover for Reference 3-937 & 3-938	8" x 5" x 7/ ₁₆ " (203 mm x 127 mm x 11 mm)
3-940	Dome Cath Cover for Reference 3-937 & 3-938	8" x 5" x 7/ ₁₆ " (203 mm x 127 mm x 11 mm)
3-941	Cath Tray, Narrow, 8 x 3 x 1 1/2	8 ³⁵ / ₆₄ " x 3 ⁵ / ₆₄ " x 1 ¹ / ₂ " (217 mm x 78 mm x 38 mm)
3-942	Dome Cath Cover, Narrow for Reference 3-941	8 ²¹ / ₃₂ " x 3 ⁵ / ₁₆ " x 19/ ₃₂ " (220 mm x 84 mm x 15 mm)
3-943	Cath Tray, Long, 15 ³ / ₄ " x 2 ³ / ₄ " X 2	15 ³ / ₄ " x 2 ³ / ₄ " x 2" (400 mm x 70 mm x 51 mm)
3-944	Dome Cath Cover, Long for Reference 3-943	15 ¹³ / ₁₆ " x 3 ¹ / ₃₂ " x 29/ ₃₂ " (402 mm x 77 mm x 23 mm)

For more information or to place an order, please contact:

Integra ▪ 589 Davies Drive, York, PA 17402

866-854-8300 USA ▪ 717-840-9335 outside USA ▪ 717-840-9347 fax

integralife.com/integra-miltex



A' [http://www.dermatology.ucsf.edu/education_training/Residency Program/SURGERY/FINAL PROCEDURE MOHS Chapter MARCH 05.pdf](http://www.dermatology.ucsf.edu/education_training/Residency%20Program/SURGERY/FINAL%20PROCEDURE%20MOHS%20Chapter%20MARCH%2005.pdf)

B' [http://www.mohssurgery.org/files/conf/papers/2009 FMS Physician/1G13_50 - Yob office setup 2009.pdf](http://www.mohssurgery.org/files/conf/papers/2009%20FMS%20Physician/1G13_50%20Yob%20office%20setup%202009.pdf)

Due to the hand craftsmanship of the instruments, tolerances in product specifications are subject to minor variations in dimensions and finishes.

Integra, the Integra logo, Miltex, Padgett, Miltex Biopblade and EZ-Zyme are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries. Copyright ©2014 Integra LifeSciences Corporation. All rights reserved. 0285000-1