

# DATASHEET: MRS42

## Mobile Respiratory Cart

Technical Data	
Height:	51"
Width:	34"
Depth:	34"
Weight:	32 lbs
Ship Weight:	34 lbs
Shipping Method:	Common Carrier

### Features:

- Complete mobile respiratory station
- Holds all equipment necessary for a ventilator patient
- Oxygen concentrator fits conveniently beneath bottom shelf
- Wheeled oxygen concentrator moves with cart
- Velcro strap holds concentrator under cart while moving
- Injection molded shelves have raised lips to keep equipment safe
- Oxygen E-tank holder optional
- Handy push bar
- 30" clearance under lower shelf
- Easy-roll 5" casters (2 locking)



Replacement Parts
LC LPCSET 5 – Replacement Casters
240 SHELF – Replacement Shelf

[Link to website](#)

**INNOVATIVE PRODUCTS UNLIMITED INC.**

2120 Industrial Drive - Niles, MI 49120 - Phone: 800-833-2826 - [www.ipu.com](http://www.ipu.com)

