



## Mayo Stand

4368

1/EA

### Product Features



**WARNING:** This product can expose you to chemicals including chromium (hexavalent compounds), which are known to the State of California to cause cancer and birth defects or other reproductive harm. For more information go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

Features a high grade stainless steel tray (18/8, 304)

Base minimizes interference with other equipment

Secure height adjustments 34" - 53"

Tray size: 12-1/2" x 19"

Base width 20", depth: 19-1/2"

Weight Capacity Standard: 15 lbs

Replacement tray available (4264)

Depth at base: 19-1/4"

Inside Base Area: 10-1/2" at widest end and 5-1/2" at narrow end

### Attributes

<b>Brand Name</b>	Tech-Med
<b>Sterile</b>	No
<b>Weight</b>	15.52054 lbs

### Packaging

<b>Packaging</b>	1/EA
<b>Case Packaging</b>	1/EA
<b>Total Case Quantity</b>	1
<b>Case Packaging Barcode</b>	(01)20665973022301