# **DD-770 & DD-770**DigiDop II Display Doppler Systems









# Newman ) Medical

Employee Owned & Operated

800-267-5549 WWW.NEWMAN-MEDICAL.COM

## **You Deserve a Better Doppler**

Newman Medical offers a full line of Doppler systems with the performance and durability you expect to meet your application - obstetrical or vascular. The DigiDop 770 or 770R Display Digital Doppler have unique features not found in other Doppler systems such as Digital Sound Clarity™. The DigiDop 701 adds the Smart Recharge System™ that reduces overall costs and reduces waste. Each DigiDop is designed with the sensitivity you need while also designed to be tough and last through years of everyday office use. All our products are Made in the USA in Colorado. Each Newman Medical Doppler and probe must pass a series of comprehensive tests to assure quality and the durability of each product.

#### Probe choices include:

- **→**D2 2MHz obstetrical probe later term
- **→**D3 3MHz obstetrical probe general purpose
- **→**D5 5MHz vascular probe deeper vessels
- **→** D8 8MHz vascular probe general purpose
- **D2W** 2MHz Waterproof obstetrical water applications
- **→**D3W 3MHz Waterproof obstetrical water applications
- **→** DPPG Audio PPG vascular PPG digit probe/sensor

### **Features:**

- Better Sound Digital Sound Clarity™ produces up to 6x signal of other Dopplers
- → Better Warranty 5 year honest warranty all parts & labor. Includes rechargeable batteries on rechargeable units.
- → Better Support & Service Knowledgeable Staff & Guaranteed Service
- **→ Smart Recharge System** (Model 701) A significant reduction in lifetime battery costs and also a reduction in waste
- → 'Unbreakable' magnetic probe holder
- → All probes are interchangeable
- → User replaceable cord eliminates service downtime
- → Powered by 3 AA batteries extended battery life of 900 mins
- → Dimensions (h x w x l): 150 x 65 x 35 mm (6.0 x 2.5 x 1.25 in)
- → Weight: 340 grams (12 ounces)
- → Shipping Information: 2 pounds in 12"x9"x3" box