

NEOMEDIC

POWERPRESS

Pneumatic Sequential Compressor



**Neomedic is a manufacturer of
air compression & vacuum therapy products**

We proudly produce air compression and vacuum products. Our products include pneumatic sequential compression circulator, DVT (deep vein thrombosis) prophylaxis, Arterial flow enhancement system, A-Venous Foot Pump, Air compressor, Aspirator (800cc ~ 3000cc), Vacuum cupping massage therapy, Microdermabrasion etc. As an established company, we offer you Air Compression and Vacuum Therapy Solutions!

Powerpress Unit



- Pressure Gauge
0~160 mmHg
- Pressure
Controller
- On / Off
Switch

Item# P1000 / Pneumatic Sequential Compressor
(HCPCS: E0651)

P 1000 Series

Cycle Time:

- Inflation Time 36 sec / Deflation Time 24 sec / Each Chamber 9 sec
- Size: 12" (L) x 8" (W) x 5.5" (H)
- Weight: 5 lbs

Indications:

- Primary Lymphedema
- Secondary lymphedema
- Venous Insufficiency
- Venous stasis ulcers
- Dysfunction of the muscle pump



Item# P1000-G / Gradient Sequential Compressor



Item# P1000-CG /
Calibrated Gradient Sequential Compressor (HCPCS: E0652)
[Each 4 chamber pressure controller]

Item# P1000-W / Intermittent Pneumatic Compressor
(HCPCS: E0650)

Aspirator



Item# P4000

- On / Off Switch

- w/800cc Canister, Hydrophobic Filter, 6' Connection Tubing
- Operating Vacuum Range: 0-26 inHg (0-660 mmHg)
- Size: 12"(L) x 8"(W) x 5.5"(H)
- Weight: 5 lbs




Item# P4000-H

- Available 800cc, 1200cc, 2000c, 3000cc Canister w/ Hydrophobic Filter, 6' Connection Tubing
- Size: 18.5"(L) x 35"(W) x 28"(H)
- Weight: 14 lbs

4 Chamber Sequential Garments

- * Effective lymphatic drainage, blood circulation without backflow or gaps
- * Seamless aircell (Between each chamber)
- * Latex Free * Non Stretch Fabric



DVT Prophylaxis Therapy System



- Power On Indicator Light
- Low Pressure Indicator Light



Side View

Item# P2000 (Foot & Calf/ Thigh)

Item# P2000-F (Foot Only)

Item# P2000-C (Calf/ Thigh Only)

- **Indication:** Deep Vein Thrombosis (DVT)
- **Cycle Time:** Inflation Time 12 sec / Deflation Time 48 sec
- **Calf/ Thigh Pressure:** 40 mmHg
- **Foot Pressure:** 80 mmHg
- **Interval Time:** 15 sec between Left & Right

- Visual and Audible Alarm
- Pump Vibrates quietly
- Can be used with just one garment
- Hook and Loop fastener: easy to fit snugly and adjust to different leg sizes
- Snap-lock (self-sealing) connectors garments to the pump
- Breathable garment prevents sweating
- Size: 12”(L) x 8”(W) x 5.5”(H) • Weight: 5 lbs

Arterial Flow Enhancement Therapy

A-Venous Foot Pump



Item# P3000

[Gradient Sequential Compressor]

- Cycle Time: Inflation Time 5 sec / Deflation Time 15 sec
Inflation each chamber 2.5 sec
- Foot Pressure: 120 mmHg • Calf Pressure: 108 mm Hg



Item# P3000-F

- Cycle Time: Inflation Time 5 sec / Deflation Time 15 sec
- Foot Pressure: 120 mmHg

Item# P3000 & Item# P3000-F

- Rapid & High Inflation/ Deflation
- Visual and Audible Alarm
- Pump Vibrates quietly
- Can be used with just one garment
- Hook and Loop fastener: easy to fit snugly and adjust to different leg sizes
- Snap-lock (self-sealing) connectors garments to the pump
- Size: 12”(L) x 8”(W) x 5.5”(H)
- Weight: 5 lbs

Indication:

- Ulcers
- Vascular Insufficiency
- Peripheral Arterial Disease(PAD)
- Limb Salvage
- Intermittent Claudication
- Critical Limb Ischemia(CLI)

Elec Auto Cupping Therapy



- Item# P5000-P**
- Simple to Use
 - Non Adjustable Intensity



- Item# P5000**
- Adjustable Suction Intensity
 - Intensity Lock System
 - Size: 12" (L) x 8"(W) x 5.5" (H)



- Item# P5000-H**
- Vacuum Gauge • Power on Indicator
 - Adjustable Suction Intensity
 - Hydrophobic Filter
 - Size: 18.5" (L) x 35"(W) x 28" (H)
 - Weight: 14 lbs



Premium Cupping Cups included

- Boilable • Total 17 Cups • 10 Pressure Pin Included
- 7 different size and Type Cups
- 1 Large Massage Cup and 1 Oval size Cup Included

Rolling Tray



Item# P7000-A



Item# P7000-B

- 13" x 18" Stain-resistant wide tray
- Five star base for total stability and durability
- Adjustable height: 30" to 42"
- Easy to clean



9601 Owensmouth Ave. #8
Chatsworth, CA. 91311

Tel: 818-998-1023 / Fax:818-998-0277
Toll Free: 866-990-1168

info@air1000.com
www.air1000.com