

Nitro Aluminum Rollator, Tall Height, 10" Casters

HCPCS: E0143 E0156-1,E0143+E0156



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Tall height ideal for users 6'0" and above
- Attractive, Euro-style design
- Brake cable inside frame for added safety
- Handle height easily adjusts with unique push button (Figure A)
- Back support height easily adjusts with tool-free thumb screw (Figure B)
- Caster fork design enhances turning radius
- Large 10" front casters allow optimal steering and rolling comfort
- Easily folds with one hand to ultra-compact size for storage (Figure C)
- Plastic clip holds the Nitro closed when frame is folded (Figure D)
- Cross-brace design allows for side-to-side folding and added stability
- Seat is durable and comfortable (Figure E)
- Handsome and removable zippered storage bag with unique attachment to keep bag securely in place when Nitro is open or folded (Figure F)
- Lightweight, aluminum frame

PRODUCTS

Description

| | |
|--------------|--------------------------------------|
| RTL10266BK-T | Nitro Rollator,Tall, Black frame,1EA |
| RTL10266-T | Rollator NITRO size L Color: red |

SPECIFICATIONS

Description

| | |
|------------------------|--|
| Carton Shipping Weight | 30 lbs |
| Handle Height | 36" - 41" |
| Seat Dimensions | 18" (W) x 23.5" (H) x 8.5" (D) |
| Seat Height | 23.6" |
| Unit Dimensions | 23" (W) x 27.75" (L) |
| Warranty | Frame: Limited Lifetime Brake and Handle Housing: 5 Year Limited |
| Weight | 17.5 lbs |
| Weight Capacity | 300 lbs |

ACCESSORIES

Description

| | |
|----------|----------------------------|
| 10266-CH | Nitro Cup Holder, 1/ea, EA |
|----------|----------------------------|

PARTS

Description

| | |
|--------------|--|
| 1026627-TALL | CrossBraceAssemblyforNitroTall, 1/ea, EA |
| 1026634 | CrossFrameEndCapsToteBagW/XFrameAtt,1/ea, EA |
| 1026601-RNR | HandBrakeAssembly,Right,10266/1,ea, EA |
| 1026601-GR | HandBrakeAssmby,Right,10266BK/10266BK-T, EA |

| PARTS | Description |
|------------------|--|
| 1026633-R | ToteBag,CrossFrameAttachment,10266,1/ea, EA |
| 1026633-BK | ToteBag,CrossFrameAttachment,10266BK,1ea, EA |
| 1026608-TALL | SeatUpholW/AttachingHardwareForNitroTall, EA |
| 1026609-FSW-TALL | SeatLiftingStrapW/AttachHrdwre,10266WT-, EA |
| 1026609-FSR-TALL | SeatLiftingStrapW/AttachHrdwre,10266-T, EA |
| 1026609-FSG-TALL | SeatLiftingStrapW/AttachHrdwre,10266BK-, EA |
| 1026601-WR | HandBrakeAssmby,Right,10266WT, 1/ea, EA |
| 1026632 | Brake Cable Only, Nitro, 1/ea, EA |
| 1026601-GR-T | HandBrakeAssmby,Right,NitroTall,Black,1, EA |
| 1026601-G-T | HandBrakeAssmby,Left,NitroTall,Black,1e, EA |
| 1026601-WR-T | HandBrakeAssmby,Right,NitroTall,White,1, EA |
| 1026601-RNR-T | HandBrakeAssmby,Right,NitroTall,Red,1e, EA |
| 1026601-RN-T | HandBrakeAssmby,Left,NitroTall,Red,1e, EA |
| 1026632-T | Brake Cable & Sleeve, Nitro Tall, 1/ea, EA |
| 1026633-WT | ToteBag,CrossFrameAttachment,10266WT,1ea, EA |
| 1026646 | CrossFrameBolt,Nitro/NitroHemiHD/NitroH, EA |