

Nitro HD Rollator

HCPCS:




FIGURE A



FIGURE B



FIGURE C



FIGURE D




FIGURE E



FIGURE F

- Attractive, Euro-style design
- Brake cable inside frame for added safety (A)
- Handle height easily adjusts with unique push button
- Unique suspension system for support and comfort (B)
- Back support height easily adjusts with tool-free thumb screw
- Caster fork design enhances turning radius (C)
- Large 10" front casters allow optimal steering and rolling comfort (C)
- Easily folds with one hand to ultra-compact size for storage (D)
- Under seat locking mechanism holds the Nitro closed when frame is folded (E)
- Cross-brace design allows for side-to-side folding and added stability
- Seat is durable and comfortable (F)
- Handsome and removable zippered storage bag with unique attachment to keep bag securely in place when Nitro is open or folded
- Lightweight, aluminum frame

PRODUCTS

Description

| | |
|------------|--------------------------------|
| 10266HD-BK | HD Nitro Rollator, Black, 1/ea |
|------------|--------------------------------|

SPECIFICATIONS

Description

| | |
|-----------------|---|
| Handle Height | 36" - 41" |
| Seat Dimensions | 23" (W) x 10" (D) |
| Seat Height | 23.5" |
| Unit Dimensions | 28" (L) x 28.25" (W) |
| Warranty | Frame: Limited Lifetime Brake Handle and Housing: 5 Years |
| Weight | 23 lbs |
| Weight Capacity | 450 lbs |

ACCESSORIES

Description

| | |
|-----------------|--|
| 10266-BRAKE ENH | Nitro Slow Down Brake, ENH, 1PR, 25/CS, EA |
|-----------------|--|

PARTS

Description

| | |
|-------------|---|
| 1026646 | CrossFrameBolt,Nitro/NitroHemiHD/NitroH, EA |
| 10266HD-33 | Cross Frame, Nitro HD, 1/ea, EA |
| 10266HD01-R | Hand Brake Assmby, Right, Nitro HD, 1ea, EA |
| 1026617-AN | FrontWheelW/Bearings,NitroW/NewFork,1/e, EA |
| 1026611 | ToteBagMountSetW/Screws for Nitro, 1/pr, PR |
| 1026641 | WheelBrakeLeftW/Hrdwr,Plastic,Nitro,1ea, EA |

| PARTS | Description |
|-------------|--|
| 1026625-AN | Wheel&ForkAssbly,Left,Alum,NitroNew,1/ea, EA |
| 1026624-AN | Wheel&ForkAssbly,Rght,Alum,NitroNew,1/ea, EA |
| 10266HD01-L | Hand Brake Assmbly, Left, Nitro HD, 1/ea, EA |
| 1026642 | WheelBrakeRightW/Hrdwr,Plastic,Nitro,1e, EA |
| 1026644 | Folding Lock Mechanism, Red, Nitro,1/ea, EA |
| 1026618-AN | FrontWheelBearings,NitroW/NewFork, 1/pr, PR |
| 1026614-AN | Hub Cap for Aluminum NitroAlum.Fork,1/ea, EA |
| 1026622-AN | Fork, Left, Aluminum, Nitro New, 1/ea, EA |
| 1026643 | Cross Brace Screw & Insert, Nitro, 1/ea, EA |
| 1026621-AN | Fork, Right, Aluminum, Nitro New, 1/ea, EA |
| 1026645 | Front Axle & Washer,NitroW/NewFork, 1/e, EA |
| 892 | Cane Holder, Nitro Rollator, 1/ea, EA |
| 1026637 | Step Tube, Right, Nitro, 1/ea, EA |
| 10266HD31 | Shock Absorber, Nitro HD, 1/ea, EA |