

## Otter Bathing System

HCPCS: E1399-1



FIGURE A




FIGURE B




FIGURE C

- Plastic frame with a choice of standard (Figure B) or soft (Figure C) removable and washable fabric
- Standard Fabric: Vinyl covered nylon dries easily (Figure B)
- Soft Fabric: Polyester knit fabric includes a trunk strap to assist in lateral support
- Features additional padding between the fabric and Otter frame (Figure C)
- Angle-adjustable seat and back have five angle adjustments at 0°, 22.5°, 45°, 67.5° and 90°
- Seat and back Uni-Bars provide one-handed angle adjustments (Figure A)
- Adjustments can be made with user seated in chair
- Height- and width-adjustable head cushions offer support and comfort
- They can also be used as lateral supports
- Removable positioning straps can be placed anywhere along the length of the frame
- Leg straps control abduction and adduction
- Adjustable slip-resistant legs raise chair up to 7"
- Folds flat for easy storage

### PRODUCTS

#### Description

OT 1000	Otter bath chair,small,1/cs
OT 2000	Otter bath chair,medium,1/cs
OT 3000	Otter bath chair,large,1/cs
OT 4000	Otter Bath Chair, Standard Fabric, XL
OT-1000SF-BB	Otter Bath Chair, Soft Fabric, S, Blue
OT-1000SF-LL	Otter Bath Chair, Soft Fabric, S, Lav
OT-2000SF-BB	Otter Bath Chair, Soft Fabric, M, Blue
OT-2000SF-LL	Otter Bath Chair, Soft Fabric, M, Lav
OT-3000SF-BB	Otter Bath Chair, Soft Fabric, L, Blue
OT-3000SF-LL	Otter Bath Chair, Soft Fabric, L, Lav
OT-4000SF-BB	Otter Bath Chair, Soft Fabric, XL, Blue
OT-4000SF-LL	Otter Bath Chair, Soft Fabric, XL, Laven

### SPECIFICATIONS

#### Description

Overall Outside Width	17"
Seat Dimensions	14.5" (W) x 2"-7" (H)
Warranty	Limited Lifetime

### ACCESSORIES

#### Description

OT 8010	Otter Tub Stand, 1/cs, CS
OT 8000	Lateral Support/Headrest,Otters,1ea, EA

**ACCESSORIES****Description**

OT 8020	Otter ShowerStand,wheeled,1/cs, CS
---------	------------------------------------

**PARTS****Description**

OT 1002	Fabric,Soft,w/allStraps,OT1000,1/bx, BX
---------	---

OT 2002	Fabric,Soft,w/allStraps,OT2000,1/bx, BX
---------	---

OT 3002	Fabric,Soft,w/allStraps,OT3000, BX
---------	------------------------------------

OT 0020	Leg Strap, OT1000/2000, 1/ea, EA
---------	----------------------------------

OT 3001	Fabric, OT3000, replcmnt, 1/ea, EA
---------	------------------------------------

OT 3014	Seat Frame, OT 3000,1/ea, EA
---------	------------------------------

OT 1BT-L	BackTrigger,Lft,OT1000,1/ea, EA
----------	---------------------------------

OT 1BT-R	Back Trigger,Rght, OT1000,1/ea, EA
----------	------------------------------------

OT 2BT-L	BackTrigger,Lft,OT2000,1/ea, EA
----------	---------------------------------

OT 2BT-R	Back Trigger, Rght,OT2000,1/ea, EA
----------	------------------------------------

OT 3BT-L	BackTrigger,Lft,OT3000,1/ea, EA
----------	---------------------------------

OT 0010	strap,positioning,OT1000/2000,1/ea, EA
---------	--

OT 3BT-R	Back Trigger,Rght, OT3000,1/ea, EA
----------	------------------------------------

OT 0011	strap,positioning,OT3000,1/ea, EA
---------	-----------------------------------

OT 1ST-L	SeatTrigger,Left,OT1000 & OT2000,1/ea, EA
----------	---

OT 3ST-L	SeatTrigger,Lft,OT3000,1/ea, EA
----------	---------------------------------

OT 1ST-R	Seat Trigger,Rght, OT1000 & OT2000,1/ea, EA
----------	---

OT 3ST-R	Seat Trigger, Rght,OT3000,1/ea, EA
----------	------------------------------------

OT 0013	Otter Unibar Mntng Bracket, 1ea, EA
---------	-------------------------------------

OT 0012	Unibar w/Mounts, Otter Back, 1/ea, EA
---------	---------------------------------------