



## Hydraulic Patient Lift

Designed and created to help make patient handling safer for caregivers by helping to prevent back injuries from manual transfers, the **ProBasics Hydraulic Patient Lift** is safe, stable and affordable. With its slim design, the ProBasics Hydraulic lift is versatile enough for use with most patient transfers. The lift features include durable steel construction with a hammertone finish, a push handle, a padded swivel bar and two locking rear casters. The ProBasics Hydraulic patient lift is light-weight, assembles and disassembles easily, and will rotate patients without side-to-side sway. Supports patient weights up to 450 pounds.



Adjustable width base with smooth low friction casters with two locking rear casters.



Rotating handle on hydraulic pump and padded bar.



Offset mast with boom style provide better lift and range.

### FEATURES

- Heavy duty round tube steel construction with a chip resistant hammertone finish
- Six-point swivel bar attachment easily adapts to all sling styles and positions (can be used with 4 or 6 point sling)
- Low friction casters significantly improve rollability and contribute to caregiver safety. Rear two casters are locking
- 17.5" - 61" range allows for patient to be picked up from lying position on floor
- Pump handle can rotate from side to side for convenience of Attendant
- Offset mast and boom style provide better lift path, maximizing full travel range
- Swivel bar and pump handle are padded for comfort
- Elongated handle design serves a broad range of operator heights and offers improved handling ability
- Completely enclosed cam moves legs of adjustable base
- Base offers obstruction-free operation; ideal for cleaning
- Screw-in pin at mast base provides maximum stability and prevents removal of mast

### SPECIFICATIONS

Overall Width of Base (when open):	42" (for maximum safety)
Overall Width of Base (when closed):	22.5" (to easily fit thru door openings)
Minimum Height:	17.5"
Maximum Height:	61"
Base Height:	6.5'
Caster Size:	5"
Weight Capacity:	450 lb
Product Weight:	72 lb
Warranty:	Product lifetime (5 years) on frame, One year on all parts
HCPCS Code:	E0630

SLINGS ON NEXT PAGE ►

Item #	Description	UOM
PBPLSM001	ProBasics Hydraulic Patient Lift	each

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



## Patient Lift Slings

### MESH SLINGS

- Four point strap hookup
- Mesh slings are washer safe and can be used in wet or dry environments
- Commode Opening: 8.5" x 11"
- Available in a variety of sizes
- HCPCS Code: E0621



4 color-coded slings per corner



3-loop patient positioning handle

Item #	Description	Dimensions	Size	Weight Capacity	UOM
<b>MESH SLINGS</b>					
PBSL110	ProBasics Full Body Mesh Sling	41.5" W x 53.3" L	Medium	450 lb	each
PBSL111	ProBasics Full Body Mesh Sling	45.5" W x 58.3" L	Large	450 lb	each
PBSL114	ProBasics Full Body Mesh Sling with Commode Opening	41.5" W x 53.3" L	Medium	450 lb	each
PBSL115	ProBasics Full Body Mesh Sling with Commode Opening	45.5" W x 58.3" L	Large	450 lb	each
PBSL116	ProBasics Full Body Mesh Sling with Commode Opening	45.5" W x 63.3" L	Extra Large	450 lb	each