



# 10233 CLEVER-LITE LS

- 1. Seat- 10230SN
- 2. Rear Glide Brake- 10233RL
- 3. Walker Tip Grey- 10325  
Black- 10325B
- 4. Walker Wheel- 10109
- 5. Adjustment Button- 12002PB
- 6. Basket- 10230B
- 7. Seat Frame- 10230SF



## SPECIFICATIONS

Item #	Width (inside hand grip)	Width (inside back legs)	Overall Width	Depth (opened at base)	Seat Width	Seat Depth	Seat Height Highest	Seat Height Lowest	Handle Ht. (highest setting)	Handle Ht. (lowest setting)	Product Weight	Carton Shipping Weight	Weight Limit
10233	18.5"	21.5"	27"	24"	19"	11"	23.5"	19.5"	39"	32"	11 lbs.	10.5 lbs.	300 lbs.