


Contracture Management > Elbow Supports

Rolyan® Rolyan® Arm Restraint System. Used to limit elbow flexion. Helps prevent clients from chewing on their hands, pulling on IV tubing or touching their face. Soft, comfortable foam-padded outside with plastic insert inside. Hook & loop closures. Sold in pairs. Latex free.



Pediatric	Size	Width x Length
802201	Small	7" x 5" (18 x 13 cm)
802202	Medium	9" x 7" (23 x 18 cm)
802203	Large	10" x 8" (25.5 x 20 cm)
Adult (solid white)		
802204	Medium	10" x 11" (25.5 x 28 cm)
802205	Large	12" x 13" (30.5 x 33 cm)

Freedom Adjustable Elbow/Knee Splints. Multipurpose soft splints help restrict elbow and knee movement or can be used to help prevent joint contractures. Filled with polyester beads for a soft way to position a limb. Machine wash. The A812 series features a zipper opening to allow you to adjust the amount of beads in the splint and is grayish in color. The 550054 series feature a nylon loop at the end of the splint that may be pinned or clipped to the patient's sleeve to prevent the splint from migrating. Latex free.



Item No.	Length	Fits Circ.
55005403	7"	3"-10"
A81240	9 1/2"	3" - 6 1/2"
55005404	10"	3"-12"
A81241	12"	7"-9"
55005405	11"	3"-14"
A81242	13 1/2"	10"-12"
A81243	13 1/2"	13"-15"
55005406	13"	6"-17"
A81244	15"	16"-18"
55005407	14"	8"-21"
55005408	18"	14"-27"


Contracture Management

The Leeder Elbow Splint. Double-upright design is perfectly suited for treatment of non-fixed contractures. Indexed hinge indicates ROM improvement, and allows dynamic or static settings. Easy to apply and comfortable to wear. Carefully controls increase of ROM and stabilizes elbow extension. 15 degrees between stops with a total of 120 degrees ROM. Latex free.



Item No.	Size	Circumference	
		Bicep	Forearm
55460801	Small	8 1/2"-10"	7 1/2"-9"
55460802	Medium	10 1/2"-13"	9 1/2"-12"

Dynamic Elbow Extension Orthosis. Increases elbow extension in patients with non-fixed contractures. Patient-friendly application. Well-padded Kydex forearm and humeral cuffs. Includes optional distal, posterior elbow strap, choice of three coil spring tensions and wrench for set screws. Latex free.



Item No.	Size	Overall Length	Bicep Circ.
081530765	Small	12" (30.5cm)	9"-11" (22.9-27.9cm)
081530773	Medium	12" (30.5cm)	11"-13" (27.9-33cm)
081530781	Large	13" (33cm)	13"-16" (33-40.6cm)