PLASTIC OVERBED TABLE

Convenient, Durable and Accessible

FEATURES AND BENEFITS

- Heavy-duty, steel-welded H-style frame
- Impermeable plastic tabletop is easy to disinfect
- Four 2.5" swivel casters for maximum maneuverability
- 2 locking casters provide added stability
- Spring-loaded locking lift mechanism adjusts table height: 28.5"-41.9"
- Available in Mahogany (10457MH), Cherry (10457CH), and Grey (10457GR)
- Distributed tabletop weight capacity 50 lb.





Reorder No.	Product Name	Max. Capacity
10457CH	Plastic-Top Overbed Table - Cherry	50 lb.
10457MH	Plastic-Top Overbed Table - Mahogany	50 lb.
10457GR	Plastic-Top Overbed Table - Grey	50 lb.







dynarex[®] is your one-stop medical products brand since 1967

For more information contact your local Dynarex sales representative OR
Call us at: (888) 396-2739 • www.dynarex.com