


Pressure Redistribution Mattress

Model # NV-PPM7



Features:

- Multi-layered Pressure Redistribution Mattress
- Economical option maximizes value and comfort with 3 full layers of foam available with or without raised side rails.
- Top layer uses unique die-cutting and a true memory foam heel section to optimize zoned pressure redistribution.
- Bottom layer provides full length horizontally scored articulation cuts extending the durability and life of the mattress significantly.
- Cover is fire retardant, fluid proof, tear resistant, low shear, and easy to clean. Bottom cover is waterproof vinyl.
- Weight Capacity: 300 lbs
- 3 Layers of Foam
- Mattress Cover Material: Fluid-Resistant Nylon Top & Waterproof Vinyl Bottom
- Fire Rating: CFR 16 Part 1633 Federal Fire Code
- Warranty: 7 Year Limited Foam, 1 Year Cover

Model Numbers

Model Number	Dimensions	Description
NV-PPM7-80	80" (L) x 36" (W) x 6" (H)	Multi-Layered/Multi-zoned Foam Mattress
NV-PPM7-84	84" (L) x 36" (W) x 6" (H)	Multi-Layered/Multi-zoned Foam Mattress
NV-PPM7-80-RS	80" (L) x 36" (W) x 6" (H)	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out
NV-PPM7-84-RS	84" (L) x 36" (W) x 6" (H)	Multi-Layered/Multi-zoned Foam Mattress w/3" Elevated Perimeter and Cut-out