Stocked First Aid Kit - 25 person

First Aid Kits - 25 person

Metal case with 1 shelf

- Ideal for contractors, fleet vehicles, work sites or small companies
- Suitable for up to 25 employees
- Enameled steel case with handle
- Can be hung on wall
- · Cover can be used as shelf when opened
- Medicine & Antiseptics
 39 BZK Antiseptic towelettes
 6 First aid/burn cream packs, 0.9 gm
 6 Antibiotic ointment pack, 0.9 gm
- Bandages

100 1" x 3" Adhesive plastic bandages

- 2 Knuckle fabric bandages
- 2 Fingertip fabric bandages
- Injury Treatment
 - 15" x 9" Trauma pad
 - 8 3" x 3" Gauze dressing pads, 4 2-pks

1 40" x 40" x 56" Triangular sling/bandage, w/2 safety pins

- $1\frac{1}{2}$ " x 10 yd. First aid tape roll
- 13" Conforming gauze roll bandage
- 1 4" x 5" Instant cold compress
- References / Instruments
 - 1 First Aid Guide
 - 4 Exam quality vinyl gloves
 - 1 41/2" Scissors, nickel plated
 - 1 4" Tweezers, plastic, one time use
 - 1 Metal case
- Total 178 pieces

Plastic case with dividers

- Ideal for contractors, fleet vehicles, work sites or small companies
- Suitable for up to 25 employees
- Water-resistant plastic case with handle and can be wall mounted
- Medicine & Antiseptics

15 BZK Antiseptic towelettes

- 6 First aid/burn cream packs, 0.9 g
- 6 Antibiotic ointment pack, .9gm
- 1 Ammonia inhalant ampoules
- Bandages

100 1" x 3" Adhesive plastic bandages

- 2 Knuckle fabric bandages
- 2 Fingertip fabric bandages
- Injury Treatment
 - 15" x 9" Trauma pad
 - $8\,3"\,x\,3"$ Gauze dressing pads, $4\,2\text{-pks}$
 - 1 40" x 40" x 56" Triangular sling/bandage, w/2 safety pins
 - $1\frac{1}{2}$ " x 10 yd. First aid tape roll
 - 13" Conforming gauze roll bandage
 - 12" Conforming gauze roll bandage
 - 3 2" x 3" Non-stick pads
 - 2 Sterile eye pads
 - 1 Evewash
 - 1 4" x 5" Instant cold compress
- References / Instruments
 - 1 First Aid Guide
 - 4 Exam quality vinyl gloves
 - 1 41/2" Scissors, nickel plated
 - 14" Tweezers, plastic, one time use
 - 1 Plastic case
- Total 162 pieces

1799-25 Metal case w/1 shelf 10 ½" x 7 ½" x 2 ½" 1 ea

1799-25P

Plastic case w/dividers 9" x 8 3/8 " x 2 1/2"

1 ea

Meets OSHA recommendations and ANSI/ISEA Z308.1-2009





