

## Patient Positioners

All of our foam positioners are:

- Latex free
- Produced with 1.5 density foam for ultimate performance and comfort
- Meet the California Technical Bulletin 117 Standard
- Non-sterile
- Cases are vacuum packaged for maximum storage capacity




*Ulna Nerve Protector with Adjustable Straps*  
**P201** 15 ½" x 6" 1 pr



*Arm Cradle*  
**P202** 12 ½" x 16 ½" x 9" 1 ea



*Heel and Elbow Protector with Adjustable Straps*  
**P204** Standard 6" x 7" x 2" 1 pr  
**P205** Boot-Shape 8" x 8 ½" x 2" 1 pr



*Abduction Pillow with Adjustable Straps*  
**P301** Small 14" x 5" x 18" 1 ea  
**P302** Medium 16" x 6" x 22" 1 ea  
**P303** Large 18" x 6" x 24" 1 ea



*Body Wedge*  
**P401** 24" x 11" x 8" 1 ea



*Foot Drop Support with Adjustable Straps*  
**P501** 7" x 8 ½" x 7" 1 ea

