



ROLLATORS

WALKABOUT CONTOUR IMPERIAL HEMI

The Walkabout ConTour Imperial Hemi is the hemi solution for the bariatric market! It has the same great features as the ConTour Imperial and folds quickly and easily for storage and travel. Features include a padded seat, ergonomic hand grips, easy-to-operate hand brakes, and a contoured, padded backrest.

Features

- 20.5" seat to floor height
- Epoxy coated aluminum frame
- Locking loop brakes
- Removable pouch under seat for storage of personal items
- Adjustable handle height
- Heavy-duty 8" wheels for indoor and outdoor use on even terrain
- Weighs 17.5 lb
- Includes laminated color operating instructions
- Designed to meet the requirements of HCPCS codes: E0149 + E0156
- Maximum Weight Capacity 450 lb - EVENLY DISTRIBUTED
- Five (5) year limited warranty



		Width inside handgrips		overall Width		overall Depth		Seat Width		Seat Depth		Seat Height	
		inch	cm	inch	cm	inch	cm	inch	cm	inch	cm	inch	cm
RJ4402R	Burgundy	21.5	54.61	28.75	73.03	26	66.04	19.5	49.53	13.8	35.05	20.5	52.07

		Wheel diameter		minimum Handle Height		maximum Handle Height		Product Weight		Weight Capacity EVENLY DISTRIBUTED		UOM
		inch	cm	inch	cm	inch	cm	lb	Kg	lb	Kg	
RJ4402R	Burgundy	8	20.32	31	78.74	34.5	87.63	17.5	7.94	450	204.12	1 ea



Rollators are not intended to be propelled while seated. Do NOT use this rollator as a wheelchair or transport device.



GF Health Products, Inc. is an
intertek ISO 13485:2016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice.
The most current and complete product information can be found on our website.

www.grahamfield.com



© 2015, GF Health Products, Inc. All Rights Reserved.
Lumex is a registered trademark of GF Health Products, Inc.

graham
field

GF-15-00039RevC22