

SUPPORT SURFACES



bariatric SKIN PROTECTION CUSHIONS

The Bariatric Skin Protection Cushions provides pressure distribution and weight equalization to meet unique demands of the bariatric patient.

Cushion Features

- 4" cushion
- Medium density foam base
 Non-aqueous specialized visco-
- elastic gel core layer reduces tissue interface pressure at critical points
- Proprietary gel is a flexible solid that will not leak
- Gel is placed at three pressure zones
- Visco-elastic memory foam topper for superior comfort and pressure relief
- Positions the user with a moderate slope from front to back

Cover Features

- Fluid-resistant stretch top
- cover, reduces shear and friction
- Antimicrobial
- Flame retardant
- Removable and washable
- Non-skid base
- Strap to secure cushion firmly to the chair

Specifications		Meets	Overall Width		Overall Length		Overall Thickness		Weight Capacity evenLy distributed		Made in US∕	Warranty in months
Item	UOM	HCPCS	inch	cm	inch	cm	inch	cm	lb	kg	Σ̈́	≥.≘
20252218B	1/ea	E2604	22	55.9	18	45.7	4	10.2	650	294.8		18
20252418B	1/ea	E2604	24	61	18	45.7	4	10.2	650	294.8		18
20252619B	1/ea	E2604	26	66	19	48.3	4	10.2	650	294.8		18
20253019B	1/ea	E2604	30	76.2	19	48.3	4	10.2	650	294.8		18





graham field

"GF Health Products, Inc. is an intertek ISO 134852016 / MDSAP Certified Company

Color variations may occur due to monitor or printer settings.

+1 770.368.4700

Information contained herein is subject to change without notice. The most current and complete product information can be found on our website. www.gendroninc.com | www.grahamfield.com

