

## Pure Expressions Dual Channel Electric Breast Pump

HCPCS: E0603-1,E0603



FIGURE A



FIGURE B



FIGURE C



FIGURE D

- LCD Display provides status of programs and functions
- Control suction and speed to suit your comfort level with a total of 28 pumping state options
- 2 Phase Expression Stimulation and Expression Mode
- Easy to assemble and clean silent for discreet pumping
- Passive prevention of milk backflow
- Single or Double Pumping Options
- BPA and Latex Free
- Two 4oz Bottle Feeding Adapter Kits and Storage Containers with Lids
- Piston Pump System imitates baby's sucking rhythm
- Includes A/C Power Supply & External battery pack 6 AA alkaline batteries (not included)
- One Year Limited Warranty

PRODUCTS	Description
RTLBP0200	Double Electric Breast Pump,1/EA
SPECIFICATIONS	Description
Batteries (pair)	6 AA Alkaline
Contains Latex	No
Product Dimensions	7" (W) x 4.75" (L) x 3.5" (H)
Product Weight	.93 lbs.
Retail Box Dimensions	10.35" (W) x 6.02" (L) x 8.66" (H)
Retail Box Weight	3 lbs.
Warranty	One Year Limited
ACCESSORIES	Description
MQ9409	BtlAdptrKitWithNipple&Cap,RTLBP0200,1/p, PKG
BP001	Breast Pump Carry Bag,1/ea, EA
BP002	Insulated Cooler Bag, 1/ea, EA
MQ9408	POWER AC ADAPTER FOR RTLBP0200, CS
MQ9130C	CAR CHARGER FOR RTLBP0200, EA
MQ9407	TUBING SET FOR RTLBP0200, EA
MQ9405	MEMBRANE FOR RTLBP0200, PR
MQ9403	SILICONE INSERT FOR RTLBP0200, PR
MQ9406	4oz STORAGE BOTTLE FOR RTLBP0200, EA
MQ9401	24MM BREAST PUMP SET FOR RTLBP0200, CS
MQ9404	VALVE SET FOR RTLBP0200, PR

ACCESSORIES

Description

MQ9402

27MM BREAST PUMP SET FOR RTLBP0200, CS