



HOT PACKS

- Contain bentonite clay filling which helps retain and evenly distribute heat
- Heavy duty nylon thread keeps them durable while 100% pure cotton keeps them soft & comfortable. Terry covers optional.



THERMALATORS

- Stainless steel units which store and heat Thermal Packs at a consistent temperature.
- Mobile and Stationary Units available in: 4 packs, 6 packs, 8 packs, and 12 packs.



SPLINT PANS

- Design with durable stainless along with simplicity and eco-friendly functionality.
- Include an overheating safety sensor, digital thermometer and ball valve drain
- Units maintain precise temperatures and give clinicians full control of splint pan settings



BOTTLE WARMERS

- Units keep ultrasonic gels and creams just above body temperature
- The –SBW option requires zero energy by simply affixing to the front of Thermalator units and absorbing excess radiated heat to warm up creams and lotions.



ADDITIONAL PRODUCTS

- Carts
- Hydro-Chlo® whirlpool antiseptic
- Whirlpool accessories
- Paraffin wax
- Bottle fillers

F-SERIES / HUBBARD TANKS / FULL BODY IMMERSION TANKS



The 425-gallon Hubbard Tank is most commonly used to treat the entire body simultaneously—arms, legs, and trunk of the body. But unlike regular hydrotherapy whirlpools, Hubbard tanks' figure-eight shape gives medical staff easier access to patients. This tank comes standard with twin turbine assemblies, and is capable of facilitating hydrotherapy as well as whole-body access.

Buy direct from
Whitehall Manufacturing:



www.rehabilitation-products.com

Find out more on
product details and offerings:



www.whitehallmfg.com



LOOKING TOWARDS OUR
75TH YEAR

CONTACT US:

Phone Number:
Tel: (800) 782-7706 or
(626) 968-6681

Fax: (626) 855-4862
Email: info@whitehallmfg.com

L1001456
1/20

Street Address:
15125 Proctor Avenue
City of Industry, CA 91746
USA

Mailing Address:
P.O. Box 3527
City of Industry, CA 91744
USA



REHABILITATION PRODUCT

SOLUTIONS FOR ALL YOUR REHABILITATION NEEDS:

- HEAT THERAPY
- COLD THERAPY
- SPLINT THERAPY
- HYDROTHERAPY
- BURN TANKS
- ADDITIONAL PRODUCTS

- All stainless steel construction, which increases longevity of the product and helps in infection prevention
- Anti-slip material adhered to tread
- Several models available
- Custom sizes available
- Bariatric capabilities

TRAINING STAIRS



- All stainless steel construction
- Both stationary & mobile whirlpools available
- Reversing turbines (for mobile models) are used to pump used water out of the tank and into an available sink or drain
- Customization available – powder coating & logos

HYDROTHERAPY



- These valves are designed to control the mixing of hot and cold water and deliver it at a constant, preset temperature
- The valve is equipped with a union end, angle stop, check valves with built in strainers and vacuum breaker
- Failure of the cold water supply will cause the hot water supply to close, eliminating any possibility of scalding

MIXING VALVES



- Larger GF-12 unit comes with stainless steel shelves for 12 standard Glacier Packs and a molded compartment for Ice Cups. Oversized packs are easily stored on the bottom shelf.
- Smaller LF-4 unit comes standard with specially designed shelves that store four included Glacier Packs, and a molded compartment for Ice Cups
- All models come with a thermometer

FREEZERS



- Constructed with a 13 mil PVC cover, our packs contain flexible filler that never freezes
- Each reusable pack is durable and provides up to 30 minutes of cold therapy treatment

COLD PACKS



Quality, Durability and Function



- 1 Splint Pan with stand, UCH-2 with the SP-1600-D
- 2 Cold Therapy with stand, GF-12 with GFS-1
- 3 Cold Pack,
- 4 Bottle Warmer, EBW-2
- 5 T-4-S
- 6 Thermalator with stand, T-12-M-LT with STR-1
- 7 Dry rack with large hook, TDR-6
- 8 Hot Pack,
- 9 Mixing Valve, MXG-15-OTG
- 10 Whirlpool Hydrolift, HYL-15-M
- 11 Whirlpool, Stationary, S-110-S
- 12 Washout hose, MXWH, optional accessory
- 13 Training Stairs, TS-200



CUSTOMIZE with your colors and logo



Shown in Enviro-Glazed® Powder-coated White and custom colors

REHABILITATION PRODUCT

Providing all your facilities needs.



Find out more on product details and offerings:
www.whitehallmfg.com
 Buy direct from Whitehall manufacturing:
www.rehabilitation-products.com

HydroChlor® Whirlpool Antiseptic
 Kills:
 Reproductive bacteria : 0.2%
 one(1) hour Tubercle bacillus : 5%
 two(2) hours Bacillus : More than ten(10) hours