

TheSharpChoice[®]

NMSDC Certified Diversity Supplier

RELI[®] PENCIL POINT Spinal Needles

Reli[®] Pencil Point Spinal Needles are the value alternative to leading brands and offer superior quality in a broad range of sizes.




- Needle gauges available: 20G to 27G
- Needle lengths available: 2.5" to 5"



- Large, clear hubs for enhanced tactile feel and visualization of CSF
- Improved patient comfort and recovery with atraumatic point
- Feature side-port dispensing
- Greatly reduce loss of CSF
- Variety of sizes for maximum flexibility
- Available in sterile and non-sterile packaging


BRAND (Company)	Large Hubs	Clear Hubs	Color-Coded Styles
Reli [®] Pencil Point Spinal Needle (MYCO Medical)	✓	✓	✓
Pencan and Sprotte (B. Braun)		✓	✓
Whitacre (BD)			✓
Whitacre and Sprotte (Kimberly Clark)		✓	✓



Large, clear hub



Precise side port



Broad range of sizes

GAUGE	LENGTH	STERILE PACKAGING		NON-STERILE PACKAGING	
		MYCO Code	Unit of Measure	MYCO Code	Unit of Measure
20G	3.5"	N/A	N/A	PP20G350	100 ea/bx
22G	2.5"	PP22G251	25 ea/bx; 4 bx/cs	PP22G250	100 ea/bx
22G	3.5"	PP22G351	25 ea/bx; 4 bx/cs	PP22G350	100 ea/bx
22G	5"	PP22G501	25 ea/bx; 4 bx/cs	PP22G500	100 ea/bx
24G	3.5"	PP24G351	25 ea/bx; 4 bx/cs	PP24G350	100 ea/bx
24G	4.75"	PP24G4341	25 ea/bx; 4 bx/cs	PP24G4340	100 ea/bx
25G	3.5"	PP25G351	25 ea/bx; 4 bx/cs	PP25G350	100 ea/bx
25G	4.75"	PP25G471	25 ea/bx; 4 bx/cs	PP25G470	100 ea/bx
25G	5"	PP25G501	25 ea/bx; 4 bx/cs	PP25G500	100 ea/bx
27G	3.5"	PP27G351	25 ea/bx; 4 bx/cs	PP27G350	100 ea/bx

As a leading NMSDC Diversity Certified Supplier, MYCO Medical is committed to "Bringing Reliability, Innovation and Value" to the Healthcare industry with an emphasis on "physician preference" surgical devices and disposables.

To order, learn more or receive complimentary evaluation samples:

Please call MYCO Medical at (800) 454-6926 or (919) 460-2535 or email us at info@mycomedical.com

See the complete line of MYCO innovative disposable and reusable medical devices at www.mycomedical.com



2015 Production Drive
Apex, NC 27539 USA

info@mycomedical.com
www.mycomedical.com

Some medical devices are not licensed for sale in Canada.