





Relief Pak® HotSpot® moist heat packs and covers

hot and cold



ideal for clinic or at-home use

- simple, effective method of applying 30 minutes of therapeutic moist heat
- soothe discomfort from stiff neck, backache, arthritis, sprains and bruises
- heat pack in water on stove, microwave, or hot pack heating unit
- when treatment is complete, store pack in a plastic bag in freezer
- always apply moist heat pack with a Relief Pak® cover or a heavy towel
- washable covers allow moist heat to pass without burning the skin;
- hook and loop fasteners on covers keep heat packs in place
- all-terry covers conform more easily to each treatment area; foam-filled covers are more economical

Relief Pak® HotSpot® moist heat packs

		each		bulk (12 each)	
packs only					
S	circular	10" diameter	11-1205 13.25	11-1205-12	142.50
T	standard	10" x 12"	11-1310 13.25	11-1310-12	142.50
U	neck	9" x 24"	11-1311 15.50	11-1311-12	167.50
V	oversize	15" x 24"	11-1312 30.00	11-1312-12	325.00
W	half size	5" x 12"	11-1313 10.25	11-1313-12	112.50
X	spine small	10" x 18"	11-1314 21.25	11-1314-12	230.00
Y	spine large	10" x 24"	11-1315 27.00	11-1315-12	290.00
Z	knee-shoulder	10" x 20"	11-1316 27.00	11-1316-12	290.00

Relief Pak® HotSpot® hot pack covers

		each		bulk (12 each)	
foam-filled terry covers					
	circular	12" diameter	11-1226 20.50	11-1226-12	222.50
	standard	20" x 24"	11-1360 20.50	11-1360-12	222.50
	neck	9" x 24"	11-1361 19.50	11-1361-12	210.00
	oversize	24.5" x 36"	11-1362 27.25	11-1362-12	295.00
	half size	12" x 15"	11-1363 17.00	11-1363-12	185.00
AB	standard with pocket		11-1364 17.75	11-1364-12	192.50
	oversize with pocket		11-1365 29.50	11-1365-12	317.50
		each		bulk (12 each)	
all-terry covers					
	circular	12" diameter	11-1227 20.50	11-1227-12	222.50
	standard	20" x 24"	11-1366 27.75	11-1366-12	300.00
	neck	9" x 24"	11-1367 21.50	11-1367-12	232.50
	oversize	24.5" x 36"	11-1368 36.50	11-1368-12	395.00
AA	half size	12" x 15"	11-1369 23.25	11-1369-12	250.00
	hand	14" x 14"	11-1376 25.25	11-1376-12	275.00
	dual hand	23" x 14"	11-1377 35.00	11-1377-12	420.00

Relief Pak® HotSpot® pack and cover sets

		each		bulk (12 each)	
pack with foam-filled terry cover sets					
	circular set	11-1308	29.50	11-1308-12	317.50
	standard set	11-1300	29.50	11-1300-12	317.50
	neck set	11-1301	33.00	11-1301-12	355.00