

## **Product features:**

- Weight capacity 250kg
- Tilt angle 0 30 degrees, infinite positioning
- Rolls over most standard toilets
- Low maintenance
- Arms removable / lockable / flip-back
- Height adjustable
- Safety mechanisms
- Footrests removable and swing away, height adjustable
- No rust aluminum frame (Grade 6061 T6)

## **Accessory items**

Refer to the accessory pages or on our website

- Stainless steel parts and hardware (Grade 304)
  - Safety belts lap & chest belts plus calf strap
- Caster Wheels 125mm x 4 locking
- Commode system

- Neck rest, fully adjustable
- Tool-free assembly





