

10 GLENSHAW STREET ORANGEBURG, NY 10962 1.888 DYNAREX | DYNAREX.COM

🖸 f 🏏 🙆 in



OFFSET BARIATRIC CANE REORDER NO. 10134 – 1 / CASE REORDER NO. 10134-6 – 6 / CASE HEAVY DUTY QUAD CANE REORDER NO. 10138 – 1 / CASE REORDER NO. 10138-2 – 2 / CASE

BARIATRIC CANES

500 LB. WEIGHT CAPACITY

- Constructed from high-quality anodized aluminum
- Contoured handle helps to reduce hand fatigue and improve balance
- 10138: Features large, four-legged base for added security
- Locking knob provides improved stability and helps to prevent rattling
- Non-skid rubber tips provide added safety
- Tool-free push-button height adjustment in 1-inch increments
- Comes in retail-friendly packaging
- Limited Lifetime Warranty
- HCPCS: E0105 (10138), E0100 (10134)



10134

10138

