## **DYNARIDE<sup>™</sup> RECLINING WHEELCHAIR** Extra Comfort and Stability

**Replacement Parts Available** 

## FEATURES AND BENEFITS

- Steel frame with powder-coated silver vein finish
- Hydraulic reclining mechanism adjusts up to 180 degrees
- Black padded, embossed vinyl upholstery
- Equipped with a padded head rest
- Detachable desk-length armrest
- Equipped with tool-free swing-away adjustable elevating leg rests
- Wheels are positioned at the back of the frame for additional stability
- Comes with anti-tippers
- Precision sealed wheel and caster bearings
- Packaged and sold as a single unit, 1 item per case



## WE CARE LIKE FAMILY<sup>™</sup>



Reorder No.	Product Name	Max. Capacity	HCPCS
10238	DynaRide Reclining Wheelchair 16"×16" Detach Desk Arm – Silver Vein, ELR	300 lb.	K0001, E1226, E0990 or K0195
10239	DynaRide Reclining Wheelchair 18"×16" Detach Desk Arm – Silver Vein, ELR	300 lb.	K0001, E1226, E0990 or K0195
10240	DynaRide Reclining Wheelchair 20"×16" Detach Desk Arm – Silver Vein, ELR	350 lb.	K0001, E1226, E0990 or K0195, E2201
10241	DynaRide Reclining Wheelchair 22"×18" Detach Desk Arm – Silver Vein, ELR	350 lb.	K0001, E1226, E0990 or K0195, E2201



## **dynarex**<sup>®</sup> is your one-stop medical products brand since 1967

For more information contact your local Dynarex sales representative OR Call us at: (888) 396-2739 • www.dynarex.com