DYNARIDE™ SERIES 1 WHEELCHAIRS

Configurable Essential Wheelchairs

Replacement Parts Available

FEATURES AND BENEFITS

- Steel frame with powder-coated silver vein finish
- Black padded nylon upholstery
- Padded armrests
- Push to lock brakes
- Equipped with tool-free adjustable swing-away footrests with heel loop
- 24" solid rubber wheels with composite hand rims
- Precision sealed wheel and caster bearings
- Packaged and sold as a single unit, 1 item per case



WE CARE LIKE FAMILY™



Reorder No.	Product Name	Max. Capacity	HCPCS
10254	DynaRide Series 1 Wheelchair - 18"×16" Seat with Fixed Full Arm, Silver Vein, FR	250 lb.	K0001, E0951
10255	DynaRide Series 1 Wheelchair - 18"×16" Seat with Detach Desk Arm, Silver Vein, FR	250 lb.	K0001, E0951



