DYNARIDE[™] SERIES 4 X-LITE WHEELCHAIRS

Extra Lightweight Wheelchairs

Replacement Parts Available

FEATURES AND BENEFITS

- · Carbon steel frame with powder-coated silver vein finish
- Black padded durable nylon upholstery (mildew and bacteria resistant)
- Flip back and removeable arms for easy transfer
- · Height-adjustable desk length arms with padded armrests
- Quick-release 24" solid rubber wheels with composite hand rims
- Triple axle provides easy transition of height to hemi level
- Push to lock brakes
- Equipped with tool-free swing-away foot rests with heel loop or adjustable elevating leg rests
- Adjustable 8" solid rubber casters
- Precision sealed wheel and caster bearings
- Packaged and sold as a single unit, 1 item per case



WE CARE LIKE FAMILY[™]



Reorder No.	Product Name	Max. Capacity	HCPCS
10274	DynaRide Series 4 X-Lite Wheelchair 18"×16" Flip Desk Arm ELR, Silver Vein	300 lb.	K0004, E0973, K0195, E0990
10275	DynaRide Series 4 X-Lite Wheelchair 20"×16" Flip Desk Arm ELR, Silver Vein	300 lb.	K0004, E0973, K0195, E0990, E2201
10276	DynaRide Series 4 X-Lite Wheelchair 18"×16" Flip Desk Arm, FR, Silver Vein	300 lb.	K0004, E0973, E0951
10277	DynaRide Series 4 X-Lite Wheelchair 20"×16" Flip Desk Arm, FR, Silver Vein	300 lb.	K0004, E0973, E0951, E2201



dynarex[®] is your one-stop medical products brand since 1967

For more information contact your local Dynarex sales representative OR Call us at: (888) 396-2739 • www.dynarex.com