

10 GLENSHAW STREET ORANGEBURG, NY 10962 1.888 DYNAREX | DYNAREX.COM

**REORDER NUMBERS:** 

10374 - FOLDING, 24"

10375 - DROP ARM, 24"

10377 - HD DROP ARM, 26"

10378 – HD DROP ARM, 30"





## BARIATRIC COMMODES

- Constructed from heavy-duty, powder-coated steel
- Height-adjustable to suit a variety of patients
- Tool-free setup allows for easy installation
- Includes commode bucket and cover
- Limited Lifetime Warranty
- HCPCS: E0168

## 10374

- Versatile, three-in-one commode: commode, raised toilet seat, or toilet safety frame
- Opens easily and folds flat
- Features convenient, removable back

## 10375 - 10378

- Easy-release drop-arm mechanism allows for safe lateral patient transfers
- Bucket slides in easily for added convenience

## WEIGHT CAPACITY

- 650 lb. (10374-10375)
- **1000 lb.** (10377-10378)





10374

10375





10377

10378

