DYNAREST™ FOAM MATTRESSES

Multi-Zone Foam Mattresses


FEATURES AND BENEFITS

- Constructed from horizontal precision-cut high-density foam
- Equipped with multiple pressure redistribution zones
- Features non-skid bottom for added security
- Removable, fluid resistant cover is easy to clean and provides low shear friction
- Anti-contamination zipper located beneath the mattress helps reduce foam core contamination
- Equipped with fire sock to meet Flame Retardant CFR 1632 and CFR 1633 standards
- Weight capacity 350 lb.

• HCPCS: E0184

dynarex WE CARE LIKE FAMILY

PERIMETER OPTIONS AVAILABLE



Reorder No.	Product Name	Quantity
10420	DynaRest Multi-Zone Foam Pressure Redistribution Mattress, 36" x 80"	1/Bx
10420P	DynaRest Multi-Zone Foam Pressure Redistribution Mattress w/Perimeter, 36" x 80"	1/Bx





For more information contact your local Dynarex sales representative OR

Call us at: (888) 396-2739 • www.dynarex.com