
## DYNAREST™ FOAM MATTRESSES

## Elite 4000 Foam Mattress


## **FEATURES AND BENEFITS**

- · Constructed from three layers of variable density foam
- Top layer features precision die cuts for improved pressure redistribution
- Memory foam heel slope zone helps reduce pressure wounds
- Stretch nylon cover is fluid-resistant, vapor permeable and provides low shear friction
- Cover features anti-contamination flap zipper
- Stretch top reduces perspiration and provides low shear friction
- Meets Flame Retardant CFR 1632 and CFR 1633 standards





WE CARE LIKE FAMILY™



HEEL SLOPE
Redistributes pressure from the vulnerable heels to the calf muscles

Reo	order No.	Product Name	Quantity
1	10429	DynaRest Elite 4000 Foam Mattress w/Visco Heel Slope, 36" x 80"	1/Bx
1	10430	DynaRest Elite 4000 Foam Mattress w/ Visco Heel Slope w/Perm, 36" x 80"	1/Bx





**dynarex**<sup>®</sup> is your one-stop medical products brand since 1967

For more information contact your local Dynarex sales representative OR

Call us at: (888) 396-2739 • www.dynarex.com