

# DYNAREST™ AIR MATTRESSES

## Airfloat Deluxe Air Mattress with Digital Pump



WE CARE LIKE FAMILY™

### FEATURES & BENEFITS

- Alternating pressure mattress replacement system
- Twenty 8-inch deep air cells made with cell-on-cell technology
- 18 laser holes provide low air loss
- Digital pressure control allows for 10, 15, 20 and 25-minute cycle times to customize patient comfort and therapy adjustments
- CPR valve allows for rapid deflation
- 3 modes: Alternating, Static, and Seat function
- Low pressure indicator light and audible alarm
- **HCPCS: E0277**



Reorder No.	Product Name	Dimensions	Cover Type	Max. Capacity
10443	DynaRest Airfloat Deluxe Air Mattress with Digital Pump	36"x80"	Stretch	400 lb.
10444	DynaRest Airfloat Deluxe with Perimeter and Digital Pump	36"x80"	Stretch	400 lb.



Find Us On:



**dynarex®** is your one-stop medical products brand since 1967

For more information contact your local Dynarex sales representative OR

Call us at: (888) 396-2739 • [www.dynarex.com](http://www.dynarex.com)